

Material Compatibility Status

Boehringer		
Equipment	Metrex Product	Compatibility Status
3800 Series/external surfaces	CaviCide™	No information
3800 Series/external surfaces	CaviWipes™	No information
3800 Series/external surfaces	Envirocide™	No information
3800 Series/external surfaces	CaviCide1™	No information
3800 Series/external surfaces	CaviWipes1™	No information
3800 Series/external surfaces	CaviWipes™ 2.0	No information
3800 Series/external surfaces	CaviWipes™ HP	No information
3800 Series/external surfaces	CaviWipes™ Bleach	Compatible
Model 3930	CaviCide™	Compatible
Model 3930	CaviWipes™	Compatible
Model 3930	Envirocide™	Compatible
Model 3930	CaviCide1™	No information
Model 3930	CaviWipes1™	No information
Model 3930	CaviWipes™ 2.0	No information
Model 3930	CaviWipes™ HP	No information
Model 3930	CaviWipes™ Bleach	No information

Equipment	Metrex Product	Compatibility Status
Model 3931	CaviCide™	Compatible
Model 3931	CaviWipes™	Compatible
Model 3931	Envirocide™	Compatible
Model 3931	CaviCide1™	No information
Model 3931	CaviWipes1™	No information
Model 3931	CaviWipes™ 2.0	No information
Model 3931	CaviWipes™ HP	No information
Model 3931	CaviWipes™ Bleach	No information