

Material Compatibility Status

Reichert		
Equipment	Metrex Product	Compatibility Status
Auto Refractor/Keratometer OptoCheK Plus	CaviCide™	No information
Auto Refractor/Keratometer OptoCheK Plus	CaviWipes™	No information
Auto Refractor/Keratometer OptoCheK Plus	Envirocide™	No information
Auto Refractor/Keratometer OptoCheK Plus	CaviCide1™	No information
Auto Refractor/Keratometer OptoCheK Plus	CaviWipes1™	No information
Auto Refractor/Keratometer OptoCheK Plus	CaviWipes™ 2.0	No information
Auto Refractor/Keratometer OptoCheK Plus	CaviWipes™ HP	Compatible
Auto Refractor/Keratometer OptoCheK Plus	CaviWipes™ Bleach	Compatible
CATS Tonometer Prism Model 12540	CaviCide™	No information
CATS Tonometer Prism Model 12540	CaviWipes™	No information
CATS Tonometer Prism Model 12540	Envirocide™	No information
CATS Tonometer Prism Model 12540	CaviCide1™	No information
CATS Tonometer Prism Model 12540	CaviWipes1™	No information

Equipment	Metrex Product	Compatibility Status
CATS Tonometer Prism Model 12540	CaviWipes™ 2.0	No information
CATS Tonometer Prism Model 12540	CaviWipes™ HP	No information
CATS Tonometer Prism Model 12540	CaviWipes™ Bleach	No information
Clear Chart 4 Model 13775/all parts	CaviCide™	No information
Clear Chart 4 Model 13775/all parts	CaviWipes™	No information
Clear Chart 4 Model 13775/all parts	Envirocide™	No information
Clear Chart 4 Model 13775/all parts	CaviCide1™	No information
Clear Chart 4 Model 13775/all parts	CaviWipes1™	No information
Clear Chart 4 Model 13775/all parts	CaviWipes™ 2.0	No information
Clear Chart 4 Model 13775/all parts	CaviWipes™ HP	Compatible
Clear Chart 4 Model 13775/all parts	CaviWipes™ Bleach	Compatible
Clear Chart 4 Model 13785/all parts	CaviCide™	No information
Clear Chart 4 Model 13785/all parts	CaviWipes™	No information
Clear Chart 4 Model 13785/all parts	Envirocide™	No information
Clear Chart 4 Model 13785/all parts	CaviCide1™	No information
Clear Chart 4 Model 13785/all parts	CaviWipes1™	No information

Equipment	Metrex Product	Compatibility Status
Clear Chart 4 Model 13785/all parts	CaviWipes™ 2.0	No information
Clear Chart 4 Model 13785/all parts	CaviWipes™ HP	Compatible
Clear Chart 4 Model 13785/all parts	CaviWipes™ Bleach	Compatible
Clear Chart 4 Model 13795/all parts	CaviCide™	No information
Clear Chart 4 Model 13795/all parts	CaviWipes™	No information
Clear Chart 4 Model 13795/all parts	Envirocide™	No information
Clear Chart 4 Model 13795/all parts	CaviCide1™	No information
Clear Chart 4 Model 13795/all parts	CaviWipes1™	No information
Clear Chart 4 Model 13795/all parts	CaviWipes™ 2.0	No information
Clear Chart 4 Model 13795/all parts	CaviWipes™ HP	Compatible
Clear Chart 4 Model 13795/all parts	CaviWipes™ Bleach	Compatible
CT210 Contact Tonometer Model 12598/Exterior	CaviCide™	No information
CT210 Contact Tonometer Model 12598/Exterior	CaviWipes™	No information
CT210 Contact Tonometer Model 12598/Exterior	Envirocide™	No information
CT210 Contact Tonometer Model 12598/Exterior	CaviCide1™	No information
CT210 Contact Tonometer Model 12598/Exterior	CaviWipes1™	No information

Equipment	Metrex Product	Compatibility Status
CT210 Contact Tonometer Model 12598/Exterior	CaviWipes™ 2.0	No information
CT210 Contact Tonometer Model 12598/Exterior	CaviWipes™ HP	Compatible
CT210 Contact Tonometer Model 12598/Exterior	CaviWipes™ Bleach	No information
CT210 Measurement Prism Model 12594	CaviCide™	No information
CT210 Measurement Prism Model 12594	CaviWipes™	No information
CT210 Measurement Prism Model 12594	Envirocide™	No information
CT210 Measurement Prism Model 12594	CaviCide1™	No information
CT210 Measurement Prism Model 12594	CaviWipes1™	No information
CT210 Measurement Prism Model 12594	CaviWipes™ 2.0	No information
CT210 Measurement Prism Model 12594	CaviWipes™ HP	Compatible
CT210 Measurement Prism Model 12594	CaviWipes™ Bleach	No information
Illuminated Phoropter Model 11676/body	CaviCide™	No information
Illuminated Phoropter Model 11676/body	CaviWipes™	No information
Illuminated Phoropter Model 11676/body	Envirocide™	No information
Illuminated Phoropter Model 11676/body	CaviCide1™	No information
Illuminated Phoropter Model 11676/body	CaviWipes1™	No information

Equipment	Metrex Product	Compatibility Status
Illuminated Phorofter Model 11676/body	CaviWipes™ 2.0	No information
Illuminated Phorofter Model 11676/body	CaviWipes™ HP	No information
Illuminated Phorofter Model 11676/body	CaviWipes™ Bleach	Compatible
Illuminated Phorofter Model 11676/forehead rest/face shield	CaviCide™	No information
Illuminated Phorofter Model 11676/forehead rest/face shield	CaviWipes™	No information
Illuminated Phorofter Model 11676/forehead rest/face shield	Envirocide™	No information
Illuminated Phorofter Model 11676/forehead rest/face shield	CaviCide1™	No information
Illuminated Phorofter Model 11676/forehead rest/face shield	CaviWipes1™	No information
Illuminated Phorofter Model 11676/forehead rest/face shield	CaviWipes™ 2.0	No information
Illuminated Phorofter Model 11676/forehead rest/face shield	CaviWipes™ HP	Compatible
Illuminated Phorofter Model 11676/forehead rest/face shield	CaviWipes™ Bleach	Compatible
Illuminated Phorofter/Lens/all models	CaviCide™	No information
Illuminated Phorofter/Lens/all models	CaviWipes™	No information

Equipment	Metrex Product	Compatibility Status
Illuminated Phoropter/Lens/all models	Envirocide™	No information
Illuminated Phoropter/Lens/all models	CaviCide1™	No information
Illuminated Phoropter/Lens/all models	CaviWipes1™	No information
Illuminated Phoropter/Lens/all models	CaviWipes™ 2.0	No information
Illuminated Phoropter/Lens/all models	CaviWipes™ HP	No information
Illuminated Phoropter/Lens/all models	CaviWipes™ Bleach	No information
iPAC Pachymeter Model 16040/measurement tip	CaviCide™	No information
iPAC Pachymeter Model 16040/measurement tip	CaviWipes™	No information
iPAC Pachymeter Model 16040/measurement tip	Envirocide™	No information
iPAC Pachymeter Model 16040/measurement tip	CaviCide1™	No information
iPAC Pachymeter Model 16040/measurement tip	CaviWipes1™	No information
iPAC Pachymeter Model 16040/measurement tip	CaviWipes™ 2.0	No information
iPAC Pachymeter Model 16040/measurement tip	CaviWipes™ HP	Compatible
iPAC Pachymeter Model 16040/measurement tip	CaviWipes™ Bleach	Compatible
Lenshek Plus/Pro Model 15180	CaviCide™	No information
Lenshek Plus/Pro Model 15180	CaviWipes™	No information

Equipment	Metrex Product	Compatibility Status
Lenschek Plus/Pro Model 15180	Envirocide™	No information
Lenschek Plus/Pro Model 15180	CaviCide1™	No information
Lenschek Plus/Pro Model 15180	CaviWipes1™	No information
Lenschek Plus/Pro Model 15180	CaviWipes™ 2.0	No information
Lenschek Plus/Pro Model 15180	CaviWipes™ HP	Compatible
Lenschek Plus/Pro Model 15180	CaviWipes™ Bleach	Compatible
Near Bright Model 16239	CaviWipes1™	No information
Near Bright Model 16239	CaviWipes™ 2.0	No information
Near Bright Model 16239	CaviWipes™ HP	Compatible
Near Bright Model 16239	CaviWipes™ Bleach	Compatible
Ocular Response Analyzer Gel 16170/air tube/display	CaviCide™	No information
Ocular Response Analyzer Gel 16170/air tube/display	CaviWipes™	No information
Ocular Response Analyzer Gel 16170/air tube/display	Envirocide™	No information
Ocular Response Analyzer Gel 16170/air tube/display	CaviCide1™	No information
Ocular Response Analyzer Gel 16170/air tube/display	CaviWipes1™	No information
Ocular Response Analyzer Gel 16170/air tube/display	CaviWipes™ 2.0	No information
Ocular Response Analyzer Gel 16170/air tube/display	CaviWipes™ HP	Compatible
Ocular Response Analyzer Gel 16170/air tube/display	CaviWipes™ Bleach	Compatible

Equipment	Metrex Product	Compatibility Status
Ocular Response Analyzer Gel 16170/Body	CaviCide™	No information
Ocular Response Analyzer Gel 16170/Body	CaviWipes™	No information
Ocular Response Analyzer Gel 16170/Body	Envirocide™	No information
Ocular Response Analyzer Gel 16170/Body	CaviCide1™	No information
Ocular Response Analyzer Gel 16170/Body	CaviWipes1™	No information
Ocular Response Analyzer Gel 16170/Body	CaviWipes™ 2.0	No information
Ocular Response Analyzer Gel 16170/Body	CaviWipes™ HP	Compatible
Ocular Response Analyzer Gel 16170/Body	CaviWipes™ Bleach	Compatible
Ocular Response Analyzer Gel 16170/forehead rest	CaviCide™	No information
Ocular Response Analyzer Gel 16170/forehead rest	CaviWipes™	No information
Ocular Response Analyzer Gel 16170/forehead rest	Envirocide™	No information
Ocular Response Analyzer Gel 16170/forehead rest	CaviCide1™	No information
Ocular Response Analyzer Gel 16170/forehead rest	CaviWipes1™	No information
Ocular Response Analyzer Gel 16170/forehead rest	CaviWipes™ 2.0	No information
Ocular Response Analyzer Gel 16170/forehead rest	CaviWipes™ HP	Compatible
Ocular Response Analyzer Gel 16170/forehead rest	CaviWipes™ Bleach	Compatible

Equipment	Metrex Product	Compatibility Status
Optochek Plus	CaviCide™	No information
Optochek Plus	CaviWipes™	No information
Optochek Plus	Envirocide™	No information
Optochek Plus	CaviCide1™	No information
Optochek Plus	CaviWipes1™	No information
Optochek Plus	CaviWipes™ 2.0	No information
Optochek Plus	CaviWipes™ HP	No information
Optochek Plus	CaviWipes™ Bleach	No information
Pachymetry iPAC	CaviCide™	No information
Pachymetry iPAC	CaviWipes™	No information
Pachymetry iPAC	Envirocide™	No information
Pachymetry iPAC	CaviCide1™	No information
Lenschek Plus/Pro Model 15185	CaviCide™	No information
Lenschek Plus/Pro Model 15185	CaviWipes™	No information
Lenschek Plus/Pro Model 15185	Envirocide™	No information
Lenschek Plus/Pro Model 15185	CaviCide1™	No information
Lenschek Plus/Pro Model 15185	CaviWipes1™	No information
Lenschek Plus/Pro Model 15185	CaviWipes™ 2.0	No information
Lenschek Plus/Pro Model 15185	CaviWipes™ HP	Compatible
Lenschek Plus/Pro Model 15185	CaviWipes™ Bleach	Compatible
LongLife Project-O-Chart Model 12804	CaviCide™	No information

Equipment	Metrex Product	Compatibility Status
LongLife Project-O-Chart Model 12804	CaviWipes™	No information
LongLife Project-O-Chart Model 12804	Envirocide™	No information
LongLife Project-O-Chart Model 12804	CaviCide1™	No information
LongLife Project-O-Chart Model 12804	CaviWipes1™	No information
LongLife Project-O-Chart Model 12804	CaviWipes™ 2.0	No information
LongLife Project-O-Chart Model 12804	CaviWipes™ HP	No information
LongLife Project-O-Chart Model 12804	CaviWipes™ Bleach	No information
Master Phorofter Model Ultramatic RX Model 11625B/body	CaviCide™	No information
Master Phorofter Model Ultramatic RX Model 11625B/body	CaviWipes™	No information
Master Phorofter Model Ultramatic RX Model 11625B/body	Envirocide™	No information
Master Phorofter Model Ultramatic RX Model 11625B/body	CaviCide1™	No information
Master Phorofter Model Ultramatic RX Model 11625W/chin rest/face shield	CaviCide™	No information
Master Phorofter Model Ultramatic RX Model 11625W/chin rest/face shield	CaviWipes™	No information

Equipment	Metrex Product	Compatibility Status
Master Phorofter Model Ultramatic RX Model 11625W/chin rest/face shield	Envirocide™	No information
Master Phorofter Model Ultramatic RX Model 11625W/chin rest/face shield	CaviCide1™	No information
Master Phorofter Model Ultramatic RX Model 11625W/chin rest/face shield	CaviWipes1™	No information
Master Phorofter Model Ultramatic RX Model 11625W/chin rest/face shield	CaviWipes™ 2.0	No information
Master Phorofter Model Ultramatic RX Model 11625W/chin rest/face shield	CaviWipes™ HP	Compatible
Master Phorofter Model Ultramatic RX Model 11625W/chin rest/face shield	CaviWipes™ Bleach	Compatible
Master Phorofter Model Ultramatic RX Model 11635B/forehead rest/face shield	CaviCide™	No information
Master Phorofter Model Ultramatic RX Model 11635B/forehead rest/face shield	CaviWipes™	No information
Master Phorofter Model Ultramatic RX Model 11635B/forehead rest/face shield	Envirocide™	No information
Master Phorofter Model Ultramatic RX Model 11635B/forehead rest/face shield	CaviCide1™	No information

Equipment	Metrex Product	Compatibility Status
Master Phoropter Model Ultramatic RX Model 11635B/forehead rest/face shield	CaviWipes1™	No information
Master Phoropter Model Ultramatic RX Model 11635B/forehead rest/face shield	CaviWipes™ 2.0	No information
Master Phoropter Model Ultramatic RX Model 11635B/forehead rest/face shield	CaviWipes™ HP	Compatible
Master Phoropter Model Ultramatic RX Model 11635B/forehead rest/face shield	CaviWipes™ Bleach	Compatible
Master Phoropter Model Ultramatic RX Model 11635W/body	CaviCide™	No information
Master Phoropter Model Ultramatic RX Model 11635W/body	CaviWipes™	No information
Master Phoropter Model Ultramatic RX Model 11635W/body	Envirocide™	No information
Master Phoropter Model Ultramatic RX Model 11635W/body	CaviCide1™	No information
Master Phoropter Model Ultramatic RX Model 11635W/body	CaviWipes1™	No information
Master Phoropter Model Ultramatic RX Model 11635W/body	CaviWipes™ 2.0	No information

Equipment	Metrex Product	Compatibility Status
Master Phoropter Model Ultramatic RX Model 11635W/body	CaviWipes™ HP	Compatible
Master Phoropter Model Ultramatic RX Model 11635W/body	CaviWipes™ Bleach	Compatible
Master Phoropter/Lens/all models	CaviCide™	No information
Master Phoropter/Lens/all models	CaviWipes™	No information
Master Phoropter/Lens/all models	Envirocide™	No information
Master Phoropter/Lens/all models	CaviCide1™	No information
Master Phoropter/Lens/all models	CaviWipes1™	No information
Master Phoropter/Lens/all models	CaviWipes™ 2.0	No information
Master Phoropter/Lens/all models	CaviWipes™ HP	No information
Master Phoropter/Lens/all models	CaviWipes™ Bleach	No information
ML1 Manual Lensometer Model 15110	CaviCide™	No information
ML1 Manual Lensometer Model 15110	CaviWipes™	No information
ML1 Manual Lensometer Model 15110	Envirocide™	No information
ML1 Manual Lensometer Model 15110	CaviCide1™	No information
ML1 Manual Lensometer Model 15110	CaviWipes1™	No information

Equipment	Metrex Product	Compatibility Status
ML1 Manual Lensometer Model 15110	CaviWipes™ 2.0	No information
ML1 Manual Lensometer Model 15110	CaviWipes™ HP	No information
ML1 Manual Lensometer Model 15110	CaviWipes™ Bleach	No information
Pneumatonometer Model 30	CaviCide™	No information
Pneumatonometer Model 30	CaviWipes™	No information
Pneumatonometer Model 30	Envirocide™	No information
Pneumatonometer Model 30	CaviCide1™	No information
Pneumatonometer Model 30	CaviWipes1™	No information
Pneumatonometer Model 30	CaviWipes™ 2.0	No information
Pneumatonometer Model 30	CaviWipes™ HP	Compatible
Pneumatonometer Model 30	CaviWipes™ Bleach	No information
PSL Portable Slit Lamp Model 15090	CaviCide™	No information
PSL Portable Slit Lamp Model 15090	CaviWipes™	No information
PSL Portable Slit Lamp Model 15090	Envirocide™	No information
PSL Portable Slit Lamp Model 15090	CaviCide1™	No information
PSL Portable Slit Lamp Model 15090	CaviWipes1™	No information
PSL Portable Slit Lamp Model 15090	CaviWipes™ 2.0	No information
PSL Portable Slit Lamp Model 15090	CaviWipes™ HP	No information
PSL Portable Slit Lamp Model 15090	CaviWipes™ Bleach	No information

Equipment	Metrex Product	Compatibility Status
PT 100 Tonometer Model 13850/all parts	CaviCide™	No information
PT 100 Tonometer Model 13850/all parts	CaviWipes™	No information
PT 100 Tonometer Model 13850/all parts	Envirocide™	No information
PT 100 Tonometer Model 13850/all parts	CaviCide1™	No information
Master Phorofter Model Ultramatic RX Model 11625B/body	CaviWipes1™	No information
Master Phorofter Model Ultramatic RX Model 11625B/body	CaviWipes™ 2.0	No information
Master Phorofter Model Ultramatic RX Model 11625B/body	CaviWipes™ HP	Compatible
Master Phorofter Model Ultramatic RX Model 11625B/body	CaviWipes™ Bleach	Compatible
Master Phorofter Model Ultramatic RX Model 11625B/chin rest/face shield	CaviCide™	No information
Master Phorofter Model Ultramatic RX Model 11625B/chin rest/face shield	CaviWipes™	No information
Master Phorofter Model Ultramatic RX Model 11625B/chin rest/face shield	Envirocide™	No information
Master Phorofter Model Ultramatic RX Model 11625B/chin rest/face shield	CaviCide1™	No information

Equipment	Metrex Product	Compatibility Status
Master Phoropter Model Ultramatic RX Model 11625B/chin rest/face shield	CaviWipes1™	No information
Master Phoropter Model Ultramatic RX Model 11625B/chin rest/face shield	CaviWipes™ 2.0	No information
Master Phoropter Model Ultramatic RX Model 11625B/chin rest/face shield	CaviWipes™ HP	Compatible
Master Phoropter Model Ultramatic RX Model 11625B/chin rest/face shield	CaviWipes™ Bleach	Compatible
Master Phoropter Model Ultramatic RX Model 11625W/body	CaviCide™	No information
Master Phoropter Model Ultramatic RX Model 11625W/body	CaviWipes™	No information
Master Phoropter Model Ultramatic RX Model 11625W/body	Envirocide™	No information
Master Phoropter Model Ultramatic RX Model 11625W/body	CaviCide1™	No information
Master Phoropter Model Ultramatic RX Model 11625W/body	CaviWipes1™	No information
Master Phoropter Model Ultramatic RX Model 11625W/body	CaviWipes™ 2.0	No information
Master Phoropter Model Ultramatic RX Model 11625W/body	CaviWipes™ HP	Compatible

Equipment	Metrex Product	Compatibility Status
Master Phorofter Model Ultramatic RX Model 11625W/body	CaviWipes™ Bleach	Compatible
Pachymetry iPAC	CaviWipes1™	No information
Pachymetry iPAC	CaviWipes™ 2.0	No information
Pachymetry iPAC	CaviWipes™ HP	Compatible
Pachymetry iPAC	CaviWipes™ Bleach	Compatible
PDM digital PD Meter Model 15020	CaviCide™	No information
PDM digital PD Meter Model 15020	CaviWipes™	No information
PDM digital PD Meter Model 15020	Envirocide™	No information
PDM digital PD Meter Model 15020	CaviCide1™	No information
PDM digital PD Meter Model 15020	CaviWipes1™	No information
PDM digital PD Meter Model 15020	CaviWipes™ 2.0	No information
PDM digital PD Meter Model 15020	CaviWipes™ HP	No information
PDM digital PD Meter Model 15020	CaviWipes™ Bleach	No information
Phorofter VRx Model 16241	CaviCide™	No information
Phorofter VRx Model 16241	CaviWipes™	No information
Phorofter VRx Model 16241	Envirocide™	No information
Phorofter VRx Model 16241	CaviCide1™	No information
Phorofter VRx Model 16241	CaviWipes1™	No information
Phorofter VRx Model 16241	CaviWipes™ 2.0	No information
Phorofter VRx Model 16241	CaviWipes™ HP	Compatible

Equipment	Metrex Product	Compatibility Status
Phoroptor VRx Model 16241	CaviWipes™ Bleach	Compatible
Model 15170/all parts	CaviCide™	No information
Model 15170/all parts	CaviWipes™	No information
Model 15170/all parts	Envirocide™	No information
Model 15170/all parts	CaviCide1™	No information
Model 15170/all parts	CaviWipes1™	No information
Model 15170/all parts	CaviWipes™ 2.0	No information
Model 15170/all parts	CaviWipes™ HP	No information
Model 15170/all parts	CaviWipes™ Bleach	No information
Near Bright Model 16238	CaviCide™	No information
Near Bright Model 16238	CaviWipes™	No information
Near Bright Model 16238	Envirocide™	No information
Near Bright Model 16238	CaviCide1™	No information
Near Bright Model 16238	CaviWipes1™	No information
Near Bright Model 16238	CaviWipes™ 2.0	No information
Near Bright Model 16238	CaviWipes™ HP	Compatible
Near Bright Model 16238	CaviWipes™ Bleach	Compatible
Near Bright Model 16239	CaviCide™	No information
Near Bright Model 16239	CaviWipes™	No information
Near Bright Model 16239	Envirocide™	No information
Near Bright Model 16239	CaviCide1™	No information
PT 100 Tonometer Model 13850/all parts	CaviWipes1™	No information
PT 100 Tonometer Model 13850/all parts	CaviWipes™ 2.0	No information
PT 100 Tonometer Model 13850/all parts	CaviWipes™ HP	No information

Equipment	Metrex Product	Compatibility Status
PT 100 Tonometer Model 13850/all parts	CaviWipes™ Bleach	No information
Refractor Visutron 900/ Model 16243	CaviCide™	No information
Refractor Visutron 900/ Model 16243	CaviWipes™	No information
Refractor Visutron 900/ Model 16243	Envirocide™	No information
Refractor Visutron 900/ Model 16243	CaviCide1™	No information
Refractor Visutron 900/ Model 16243	CaviWipes1™	No information
Refractor Visutron 900/ Model 16243	CaviWipes™ 2.0	No information
Refractor Visutron 900/ Model 16243	CaviWipes™ HP	Compatible
Refractor Visutron 900/ Model 16243	CaviWipes™ Bleach	Compatible
Reichert 7 Ototonometer Model 16050/air tube	CaviCide™	No information
Reichert 7 Ototonometer Model 16050/air tube	CaviWipes™	No information
Reichert 7 Ototonometer Model 16050/air tube	Envirocide™	No information
Reichert 7 Ototonometer Model 16050/air tube	CaviCide1™	No information
Reichert 7 Ototonometer Model 16050/air tube	CaviWipes1™	No information
Reichert 7 Ototonometer Model 16050/air tube	CaviWipes™ 2.0	No information
Reichert 7 Ototonometer Model 16050/air tube	CaviWipes™ HP	Compatible

Equipment	Metrex Product	Compatibility Status
Reichert 7 Ototonometer Model 16050/air tube	CaviWipes™ Bleach	Compatible
Reichert 7 Ototonometer Model 16050/forehead rest	CaviWipes1™	No information
Reichert 7 Ototonometer Model 16050/forehead rest	CaviWipes™ 2.0	No information
Reichert 7 Ototonometer Model 16050/forehead rest	CaviWipes™ HP	Compatible
Reichert 7 Ototonometer Model 16050/forehead rest	CaviWipes™ Bleach	Compatible
Sightchek Model 16290	CaviCide™	No information
Sightchek Model 16290	CaviWipes™	No information
Sightchek Model 16290	Envirocide™	No information
Sightchek Model 16290	CaviCide1™	No information
Sightchek Model 16290	CaviWipes1™	No information
Sightchek Model 16290	CaviWipes™ 2.0	No information
Sightchek Model 16290	CaviWipes™ HP	Compatible
Sightchek Model 16290	CaviWipes™ Bleach	Compatible
Tono-Pen AVIA Model 230650	CaviCide™	No information
Tono-Pen AVIA Model 230650	CaviWipes™	No information
Tono-Pen AVIA Model 230650	Envirocide™	No information
Tono-Pen AVIA Model 230650	CaviCide1™	No information
Tono-Pen AVIA Model 230650	CaviWipes1™	No information
Tono-Pen AVIA Model 230650	CaviWipes™ 2.0	No information
Tono-Pen AVIA Model 230650	CaviWipes™ HP	Compatible
Tono-Pen AVIA Model 230650	CaviWipes™ Bleach	Compatible
Tono-Pen AVIA VET Model 230650v	CaviCide™	No information

Equipment	Metrex Product	Compatibility Status
Tono-Pen AVIA VET Model 230650v	CaviWipes™	No information
Tono-Pen AVIA VET Model 230650v	Envirocide™	No information
Tono-Pen AVIA VET Model 230650v	CaviCide1™	No information
Tono-Pen AVIA VET Model 230650v	CaviWipes1™	No information
Tono-Pen AVIA VET Model 230650v	CaviWipes™ 2.0	No information
Tono-Pen AVIA VET Model 230650v	CaviWipes™ HP	Compatible
Tono-Pen AVIA VET Model 230650v	CaviWipes™ Bleach	Compatible
Tono-Pen XI Model 230635	CaviCide™	No information
Tono-Pen XI Model 230635	CaviWipes™	No information
Tono-Pen XI Model 230635	Envirocide™	No information
Tono-Pen XI Model 230635	CaviCide1™	No information
Tono-Pen XI Model 230635	CaviWipes1™	No information
Tono-Pen XI Model 230635	CaviWipes™ 2.0	No information
Tono-Pen XI Model 230635	CaviWipes™ HP	Compatible
Tono-Pen XI Model 230635	CaviWipes™ Bleach	Compatible
Reichert 7 Ototonometer Model 16050/Body	CaviCide™	No information
Reichert 7 Ototonometer Model 16050/Body	CaviWipes™	No information
Reichert 7 Ototonometer Model 16050/Body	Envirocide™	No information
Reichert 7 Ototonometer Model 16050/Body	CaviCide1™	No information

Equipment	Metrex Product	Compatibility Status
Reichert 7 Ototonometer Model 16050/Body	CaviWipes1™	No information
Reichert 7 Ototonometer Model 16050/Body	CaviWipes™ 2.0	No information
Reichert 7 Ototonometer Model 16050/Body	CaviWipes™ HP	Compatible
Reichert 7 Ototonometer Model 16050/Body	CaviWipes™ Bleach	Compatible
Reichert 7 Ototonometer Model 16050/display	CaviCide™	No information
Reichert 7 Ototonometer Model 16050/display	CaviWipes™	No information
Reichert 7 Ototonometer Model 16050/display	Envirocide™	No information
Reichert 7 Ototonometer Model 16050/display	CaviCide1™	No information
Reichert 7 Ototonometer Model 16050/display	CaviWipes1™	No information
Reichert 7 Ototonometer Model 16050/display	CaviWipes™ 2.0	No information
Reichert 7 Ototonometer Model 16050/display	CaviWipes™ HP	Compatible
Reichert 7 Ototonometer Model 16050/display	CaviWipes™ Bleach	Compatible
Reichert 7 Ototonometer Model 16050/forehead rest	CaviCide™	No information
Reichert 7 Ototonometer Model 16050/forehead rest	CaviWipes™	No information
Reichert 7 Ototonometer Model 16050/forehead rest	Envirocide™	No information
Reichert 7 Ototonometer Model 16050/forehead rest	CaviCide1™	No information