

Material Compatibility Status

VaTech		
Equipment	Metrex Product	Compatibility Status
EZsensor HD sensor & cable	CaviCide™	No information
EZsensor HD sensor & cable	CaviWipes™	No information
EZsensor HD sensor & cable	Envirocide™	No information
EZsensor HD sensor & cable	CaviCide1™	No information
EZsensor HD sensor & cable	CaviWipes1™	No information
EZsensor HD sensor & cable	CaviWipes™ 2.0	No information
EZsensor HD sensor & cable	CaviWipes™ HP	No information
EZsensor HD sensor & cable	CaviWipes™ Bleach	No information
PaX-i	CaviCide™	No information
PaX-i	CaviWipes™	Compatible
PaX-i	Envirocide™	No information
PaX-i	CaviCide1™	No information
PaX-i	CaviWipes1™	Compatible
PaX-i	CaviWipes™ 2.0	Compatible
PaX-i	CaviWipes™ HP	No information
PaX-i	CaviWipes™ Bleach	No information
PaX-i3D	CaviCide™	No information
PaX-i3D	CaviWipes™	Compatible
PaX-i3D	Envirocide™	No information
PaX-i3D	CaviCide1™	No information
PaX-i3D	CaviWipes1™	Compatible
PaX-i3D	CaviWipes™ 2.0	Compatible

Equipment	Metrex Product	Compatibility Status
PaX-i3D	CaviWipes™ HP	No information
PaX-i3D	CaviWipes™ Bleach	No information