

# Material Compatibility Status

| Zoll Medical                |                   |                      |
|-----------------------------|-------------------|----------------------|
| Equipment                   | Metrex Product    | Compatibility Status |
| 330 Multifunction Aspirator | CaviCide™         | Compatible           |
| 330 Multifunction Aspirator | CaviWipes™        | Compatible           |
| 330 Multifunction Aspirator | Envirocide™       | Compatible           |
| 330 Multifunction Aspirator | CaviCide1™        | No information       |
| 330 Multifunction Aspirator | CaviWipes1™       | No information       |
| 330 Multifunction Aspirator | CaviWipes™ 2.0    | No information       |
| 330 Multifunction Aspirator | CaviWipes™ HP     | Compatible           |
| 330 Multifunction Aspirator | CaviWipes™ Bleach | No information       |
| M Series                    | CaviWipes1™       | No information       |
| M Series                    | CaviWipes™ 2.0    | No information       |
| M Series                    | CaviWipes™ HP     | No information       |
| M Series                    | CaviWipes™ Bleach | Compatible           |
| Propaq MD                   | CaviCide™         | No information       |
| Propaq MD                   | CaviWipes™        | No information       |
| Propaq MD                   | Envirocide™       | No information       |
| Propaq MD                   | CaviCide1™        | No information       |
| Propaq MD                   | CaviWipes1™       | No information       |
| Propaq MD                   | CaviWipes™ 2.0    | No information       |
| Propaq MD                   | CaviWipes™ HP     | Compatible           |
| Propaq MD                   | CaviWipes™ Bleach | No information       |
| R Series                    | CaviCide™         | Compatible           |
| R Series                    | CaviWipes™        | Compatible           |

| <b>Equipment</b>  | <b>Metrex Product</b> | <b>Compatibility Status</b> |
|-------------------|-----------------------|-----------------------------|
| R Series          | Envirocide™           | No information              |
| R Series          | CaviCide1™            | No information              |
| R Series          | CaviWipes1™           | No information              |
| R Series          | CaviWipes™ 2.0        | No information              |
| R Series          | CaviWipes™ HP         | Compatible                  |
| R Series          | CaviWipes™ Bleach     | Compatible                  |
| Thermoguard XP    | CaviCide™             | No information              |
| Thermoguard XP    | CaviWipes™            | No information              |
| Thermoguard XP    | Envirocide™           | No information              |
| Thermoguard XP    | CaviCide1™            | No information              |
| Thermoguard XP    | CaviWipes1™           | No information              |
| Thermoguard XP    | CaviWipes™ 2.0        | No information              |
| Thermoguard XP    | CaviWipes™ HP         | No information              |
| Thermoguard XP    | CaviWipes™ Bleach     | Compatible                  |
| X Series          | CaviCide™             | No information              |
| X Series          | CaviWipes™            | No information              |
| X Series          | Envirocide™           | No information              |
| X Series          | CaviCide1™            | No information              |
| X Series          | CaviWipes1™           | No information              |
| X Series          | CaviWipes™ 2.0        | No information              |
| X Series          | CaviWipes™ HP         | Compatible                  |
| X Series          | CaviWipes™ Bleach     | No information              |
| X Series Advanced | CaviCide™             | Compatible                  |
| X Series Advanced | CaviWipes™            | Compatible                  |
| X Series Advanced | Envirocide™           | Compatible                  |
| X Series Advanced | CaviCide1™            | No information              |
| AED 3 BLS         | CaviCide™             | No information              |

| <b>Equipment</b> | <b>Metrex Product</b> | <b>Compatibility Status</b> |
|------------------|-----------------------|-----------------------------|
| AED 3 BLS        | CaviWipes™            | No information              |
| AED 3 BLS        | Envirocide™           | No information              |
| AED 3 BLS        | CaviCide1™            | No information              |
| AED 3 BLS        | CaviWipes1™           | No information              |
| AED 3 BLS        | CaviWipes™ 2.0        | No information              |
| AED 3 BLS        | CaviWipes™ HP         | No information              |
| AED 3 BLS        | CaviWipes™ Bleach     | No information              |
| AED Plus         | CaviCide™             | No information              |
| AED Plus         | CaviWipes™            | No information              |
| AED Plus         | Envirocide™           | No information              |
| AED Plus         | CaviCide1™            | No information              |
| AED Plus         | CaviWipes1™           | No information              |
| AED Plus         | CaviWipes™ 2.0        | No information              |
| AED Plus         | CaviWipes™ HP         | No information              |
| AED Plus         | CaviWipes™ Bleach     | Compatible                  |
| AED Pro          | CaviCide™             | No information              |
| AED Pro          | CaviWipes™            | No information              |
| AED Pro          | Envirocide™           | No information              |
| AED Pro          | CaviCide1™            | No information              |
| AED Pro          | CaviWipes1™           | No information              |
| AED Pro          | CaviWipes™ 2.0        | No information              |
| AED Pro          | CaviWipes™ HP         | No information              |
| AED Pro          | CaviWipes™ Bleach     | Compatible                  |
| AutoPulse        | CaviCide™             | No information              |
| AutoPulse        | CaviWipes™            | No information              |
| AutoPulse        | Envirocide™           | No information              |
| AutoPulse        | CaviCide1™            | No information              |

| <b>Equipment</b>    | <b>Metrex Product</b> | <b>Compatibility Status</b> |
|---------------------|-----------------------|-----------------------------|
| AutoPulse           | CaviWipes1™           | No information              |
| AutoPulse           | CaviWipes™ 2.0        | No information              |
| AutoPulse           | CaviWipes™ HP         | No information              |
| AutoPulse           | CaviWipes™ Bleach     | No information              |
| AutoPulse Model 100 | CaviCide™             | Compatible                  |
| AutoPulse Model 100 | CaviWipes™            | Compatible                  |
| AutoPulse Model 100 | Envirocide™           | Compatible                  |
| AutoPulse Model 100 | CaviCide1™            | Compatible                  |
| AutoPulse Model 100 | CaviWipes1™           | Compatible                  |
| AutoPulse Model 100 | CaviWipes™ 2.0        | Compatible                  |
| AutoPulse Model 100 | CaviWipes™ HP         | Compatible                  |
| AutoPulse Model 100 | CaviWipes™ Bleach     | Compatible                  |
| E Series            | CaviCide™             | No information              |
| E Series            | CaviWipes™            | No information              |
| E Series            | Envirocide™           | No information              |
| E Series            | CaviCide1™            | No information              |
| E Series            | CaviWipes1™           | No information              |
| E Series            | CaviWipes™ 2.0        | No information              |
| E Series            | CaviWipes™ HP         | No information              |
| E Series            | CaviWipes™ Bleach     | Compatible                  |
| EMV+                | CaviCide™             | Compatible                  |
| EMV+                | CaviWipes™            | Compatible                  |
| EMV+                | Envirocide™           | Compatible                  |
| EMV+                | CaviCide1™            | No information              |
| EMV+                | CaviWipes1™           | No information              |
| EMV+                | CaviWipes™ 2.0        | No information              |
| EMV+                | CaviWipes™ HP         | Compatible                  |

| <b>Equipment</b>  | <b>Metrex Product</b> | <b>Compatibility Status</b> |
|-------------------|-----------------------|-----------------------------|
| EMV+              | CaviWipes™ Bleach     | No information              |
| M Series          | CaviCide™             | No information              |
| M Series          | CaviWipes™            | No information              |
| M Series          | Envirocide™           | No information              |
| M Series          | CaviCide1™            | No information              |
| X Series Advanced | CaviWipes1™           | No information              |
| X Series Advanced | CaviWipes™ 2.0        | No information              |
| X Series Advanced | CaviWipes™ HP         | Compatible                  |
| X Series Advanced | CaviWipes™ Bleach     | Compatible                  |
| Z Vent            | CaviCide™             | Compatible                  |
| Z Vent            | CaviWipes™            | Compatible                  |
| Z Vent            | Envirocide™           | Compatible                  |
| Z Vent            | CaviCide1™            | No information              |
| Z Vent            | CaviWipes1™           | No information              |
| Z Vent            | CaviWipes™ 2.0        | No information              |
| Z Vent            | CaviWipes™ HP         | Compatible                  |
| Z Vent            | CaviWipes™ Bleach     | No information              |