

# Material Compatibility Status

| 3Shape                         |                   |                      |
|--------------------------------|-------------------|----------------------|
| Equipment                      | Metrex Product    | Compatibility Status |
| Trios 5                        | CaviCide™         | Compatible           |
| Trios 5                        | CaviWipes™        | Compatible           |
| Trios 5                        | Envirocide™       | Compatible           |
| Trios 5                        | CaviCide1™        | No information       |
| Trios 5                        | CaviWipes1™       | No information       |
| Trios 5                        | CaviWipes™ 2.0    | Compatible           |
| Trios 5                        | CaviWipes™ HP     | No information       |
| Trios 5                        | CaviWipes™ Bleach | No information       |
| Trios 3 wireless               | CaviCide™         | No information       |
| Trios 3 wireless               | CaviWipes™        | No information       |
| Trios 3 wireless               | Envirocide™       | No information       |
| Trios 3 wireless               | CaviCide1™        | No information       |
| Trios 3 wireless               | CaviWipes1™       | No information       |
| Trios 3 wireless               | CaviWipes™ 2.0    | No information       |
| Trios 3 wireless               | CaviWipes™ HP     | No information       |
| Trios 3 wireless               | CaviWipes™ Bleach | No information       |
| Trios 3 wireless Otical window | CaviCide™         | No information       |
| Trios 3 wireless Otical window | CaviWipes™        | No information       |
| Trios 3 wireless Otical window | Envirocide™       | No information       |

| <b>Equipment</b>               | <b>Metrex Product</b> | <b>Compatibility Status</b> |
|--------------------------------|-----------------------|-----------------------------|
| Trios 3 wireless Otical window | CaviCide1™            | No information              |
| Trios 3 wireless Otical window | CaviWipes1™           | No information              |
| Trios 3 wireless Otical window | CaviWipes™ 2.0        | No information              |
| Trios 3 wireless Otical window | CaviWipes™ HP         | No information              |
| Trios 3 wireless Otical window | CaviWipes™ Bleach     | No information              |