

Material Compatibility Status

BIOMERIEUX		
Equipment	Metrex Product	Compatibility Status
bioMerieux VITEK 2	CaviWipes1™	No information
bioMerieux VITEK 2	CaviWipes™ 2.0	No information
bioMerieux VITEK 2	CaviWipes™ HP	No information
bioMerieux VITEK 2	CaviWipes™ Bleach	Compatible
NucliSENS easyMAG	CaviCide™	No information
NucliSENS easyMAG	CaviWipes™	No information
NucliSENS easyMAG	Envirocide™	No information
NucliSENS easyMAG	CaviCide1™	No information
NucliSENS easyMAG	CaviWipes1™	No information
NucliSENS easyMAG	CaviWipes™ 2.0	No information
NucliSENS easyMAG	CaviWipes™ HP	No information
NucliSENS easyMAG	CaviWipes™ Bleach	Compatible
PREVI Color Gram	CaviCide™	No information
PREVI Color Gram	CaviWipes™	No information
PREVI Color Gram	Envirocide™	No information
PREVI Color Gram	CaviCide1™	No information
PREVI Color Gram	CaviWipes1™	No information
PREVI Color Gram	CaviWipes™ 2.0	No information
PREVI Color Gram	CaviWipes™ HP	No information
PREVI Color Gram	CaviWipes™ Bleach	No information
PREVI Isola	CaviCide™	Compatible
PREVI Isola	CaviWipes™	Compatible

Equipment	Metrex Product	Compatibility Status
PREVI Isola	Envirocide™	Compatible
PREVI Isola	CaviCide1™	Compatible
PREVI Isola	CaviWipes1™	Compatible
PREVI Isola	CaviWipes™ 2.0	Compatible
PREVI Isola	CaviWipes™ HP	No information
PREVI Isola	CaviWipes™ Bleach	No information
TEMPO reading station/external	CaviCide™	No information
TEMPO reading station/external	CaviWipes™	No information
TEMPO reading station/external	Envirocide™	No information
TEMPO reading station/external	CaviCide1™	No information
TEMPO reading station/external	CaviWipes1™	No information
TEMPO reading station/external	CaviWipes™ 2.0	No information
TEMPO reading station/external	CaviWipes™ HP	No information
TEMPO reading station/external	CaviWipes™ Bleach	No information
Air IDEAL 3P	CaviCide™	Compatible
Air IDEAL 3P	CaviWipes™	Compatible
Air IDEAL 3P	Envirocide™	Compatible
Air IDEAL 3P	CaviCide1™	Compatible
Air IDEAL 3P	CaviWipes1™	Compatible
Air IDEAL 3P	CaviWipes™ 2.0	Compatible
Air IDEAL 3P	CaviWipes™ HP	No information
Air IDEAL 3P	CaviWipes™ Bleach	No information

Equipment	Metrex Product	Compatibility Status
Bact/alert 3D	CaviCide™	No information
Bact/alert 3D	CaviWipes™	No information
Bact/alert 3D	Envirocide™	No information
Bact/alert 3D	CaviCide1™	No information
Bact/alert 3D	CaviWipes1™	No information
Bact/alert 3D	CaviWipes™ 2.0	No information
Bact/alert 3D	CaviWipes™ HP	No information
Bact/alert 3D	CaviWipes™ Bleach	Compatible
bioMerieux VIDAS	CaviCide™	No information
bioMerieux VIDAS	CaviWipes™	No information
bioMerieux VIDAS	Envirocide™	No information
bioMerieux VIDAS	CaviCide1™	No information
bioMerieux VIDAS	CaviWipes1™	No information
bioMerieux VIDAS	CaviWipes™ 2.0	No information
bioMerieux VIDAS	CaviWipes™ HP	No information
bioMerieux VIDAS	CaviWipes™ Bleach	Compatible
bioMerieux VITEK 2	CaviCide™	No information
bioMerieux VITEK 2	CaviWipes™	No information
bioMerieux VITEK 2	Envirocide™	No information
bioMerieux VITEK 2	CaviCide1™	No information