

Material Compatibility Status

DCI Edge		
Equipment	Metrex Product	Compatibility Status
Orthodontic chair	CaviCide™	Compatible
Orthodontic chair	CaviWipes™	Compatible
Orthodontic chair	Envirocide™	Compatible
Orthodontic chair	CaviCide1™	Compatible
Orthodontic chair	CaviWipes1™	Compatible
Orthodontic chair	CaviWipes™ 2.0	Compatible
Orthodontic chair	CaviWipes™ HP	Not recommended for use
Orthodontic chair	CaviWipes™ Bleach	Not recommended for use