

Material Compatibility Status

DentalEZ		
Equipment	Metrex Product	Compatibility Status
Plastic & coated Metal surfaces of dental chairs	CaviCide™	Compatible
Plastic & coated Metal surfaces of dental chairs	CaviWipes™	Compatible
Plastic & coated Metal surfaces of dental chairs	Envirocide™	Compatible
Plastic & coated Metal surfaces of dental chairs	CaviCide1™	No information
Plastic & coated Metal surfaces of dental chairs	CaviWipes1™	No information
Plastic & coated Metal surfaces of dental chairs	CaviWipes™ 2.0	No information
Plastic & coated Metal surfaces of dental chairs	CaviWipes™ HP	No information
Plastic & coated Metal surfaces of dental chairs	CaviWipes™ Bleach	No information
NuSmplicity Chair	CaviCide™	Compatible
NuSmplicity Chair	CaviWipes™	Compatible
NuSmplicity Chair	Envirocide™	Compatible
NuSmplicity Chair	CaviCide1™	No information
NuSmplicity Chair	CaviWipes1™	No information
NuSmplicity Chair	CaviWipes™ 2.0	No information
NuSmplicity Chair	CaviWipes™ HP	No information
NuSmplicity Chair	CaviWipes™ Bleach	Compatible

Equipment	Metrex Product	Compatibility Status
Upholstery	CaviCide™	Compatible
Upholstery	CaviWipes™	Compatible
Upholstery	Envirocide™	Compatible
Upholstery	CaviCide1™	No information
Upholstery	CaviWipes1™	No information
Upholstery	CaviWipes™ 2.0	No information
Upholstery	CaviWipes™ HP	No information
Upholstery	CaviWipes™ Bleach	Compatible