

Material Compatibility Status

Covaris		
Equipment	Metrex Product	Compatibility Status
cryoPREP Extraction system	CaviCide™	No information
cryoPREP Extraction system	CaviWipes™	No information
cryoPREP Extraction system	Envirocide™	No information
cryoPREP Extraction system	CaviCide1™	No information
cryoPREP Extraction system	CaviWipes1™	No information
cryoPREP Extraction system	CaviWipes™ 2.0	No information
cryoPREP Extraction system	CaviWipes™ HP	No information
cryoPREP Extraction system	CaviWipes™ Bleach	Compatible