

Material Compatibility Status

Gendex		
Equipment	Metrex Product	Compatibility Status
GXS-700	CaviCide™	Compatible
GXS-700	CaviWipes™	Compatible
GXS-700	Envirocide™	Compatible
GXS-700	CaviCide1™	No information
GXS-700	CaviWipes1™	No information
GXS-700	CaviWipes™ 2.0	No information
GXS-700	CaviWipes™ HP	No information
GXS-700	CaviWipes™ Bleach	No information
GXS sensor & cable	CaviCide™	Compatible
GXS sensor & cable	CaviWipes™	Compatible
GXS sensor & cable	Envirocide™	Compatible
GXS sensor & cable	CaviCide1™	No information
GXS sensor & cable	CaviWipes1™	No information
GXS sensor & cable	CaviWipes™ 2.0	No information
GXS sensor & cable	CaviWipes™ HP	No information
GXS sensor & cable	CaviWipes™ Bleach	No information