

## Material Compatibility Status

| Acteon    |                   |                      |
|-----------|-------------------|----------------------|
| Equipment | Metrex Product    | Compatibility Status |
| Mini LED  | CaviCide™         | No information       |
| Mini LED  | CaviWipes™        | No information       |
| Mini LED  | Envirocide™       | No information       |
| Mini LED  | CaviCide1™        | No information       |
| Mini LED  | CaviWipes1™       | No information       |
| Mini LED  | CaviWipes™ 2.0    | No information       |
| Mini LED  | CaviWipes™ HP     | No information       |
| Mini LED  | CaviWipes™ Bleach | No information       |
| Sopro 617 | CaviCide™         | No information       |
| Sopro 617 | CaviWipes™        | No information       |
| Sopro 617 | Envirocide™       | No information       |
| Sopro 617 | CaviCide1™        | No information       |
| Sopro 617 | CaviWipes1™       | No information       |
| Sopro 617 | CaviWipes™ 2.0    | No information       |
| Sopro 617 | CaviWipes™ HP     | No information       |
| Sopro 617 | CaviWipes™ Bleach | No information       |
| Sopro 717 | CaviCide™         | No information       |
| Sopro 717 | CaviWipes™        | No information       |
| Sopro 717 | Envirocide™       | No information       |
| Sopro 717 | CaviCide1™        | No information       |
| Sopro 717 | CaviWipes1™       | No information       |

| <b>Equipment</b>            | <b>Metrex Product</b> | <b>Compatibility Status</b> |
|-----------------------------|-----------------------|-----------------------------|
| Sopro 717                   | CaviWipes™ 2.0        | No information              |
| Sopro 717                   | CaviWipes™ HP         | No information              |
| Sopro 717                   | CaviWipes™ Bleach     | No information              |
| SoproCare/handpiece         | CaviCide™             | No information              |
| SoproCare/handpiece         | CaviWipes™            | No information              |
| SoproCare/handpiece         | Envirocide™           | No information              |
| SoproCare/handpiece         | CaviCide1™            | No information              |
| SoproCare/handpiece         | CaviWipes1™           | No information              |
| SoproCare/handpiece         | CaviWipes™ 2.0        | No information              |
| SoproCare/handpiece         | CaviWipes™ HP         | No information              |
| SoproCare/handpiece         | CaviWipes™ Bleach     | No information              |
| Sopix Series sensor & cable | CaviCide™             | Compatible                  |
| Sopix Series sensor & cable | CaviWipes™            | Compatible                  |
| Sopix Series sensor & cable | Envirocide™           | Compatible                  |
| Sopix Series sensor & cable | CaviCide1™            | No information              |
| Sopix Series sensor & cable | CaviWipes1™           | No information              |
| Sopix Series sensor & cable | CaviWipes™ 2.0        | No information              |
| Sopix Series sensor & cable | CaviWipes™ HP         | No information              |
| Sopix Series sensor & cable | CaviWipes™ Bleach     | No information              |
| Sopix Series surface        | CaviCide™             | Compatible                  |
| Sopix Series surface        | CaviWipes™            | Compatible                  |
| Sopix Series surface        | Envirocide™           | Compatible                  |
| Sopix Series surface        | CaviCide1™            | No information              |
| Sopix Series surface        | CaviWipes1™           | No information              |
| Sopix Series surface        | CaviWipes™ 2.0        | No information              |
| Sopix Series surface        | CaviWipes™ HP         | No information              |

| <b>Equipment</b>     | <b>Metrex Product</b> | <b>Compatibility Status</b> |
|----------------------|-----------------------|-----------------------------|
| Sopix Series surface | CaviWipes™ Bleach     | No information              |
| SoproCare            | CaviCide™             | No information              |
| SoproCare            | CaviWipes™            | No information              |
| SoproCare            | Envirocide™           | No information              |
| SoproCare            | CaviCide1™            | No information              |
| SoproCare            | CaviWipes1™           | No information              |
| SoproCare            | CaviWipes™ 2.0        | No information              |
| SoproCare            | CaviWipes™ HP         | No information              |
| SoproCare            | CaviWipes™ Bleach     | No information              |
| Sopros 617           | CaviCide™             | No information              |
| Sopros 617           | CaviWipes™            | No information              |
| Sopros 617           | Envirocide™           | No information              |
| Sopros 617           | CaviCide1™            | No information              |
| Sopros 617           | CaviWipes1™           | No information              |
| Sopros 617           | CaviWipes™ 2.0        | No information              |
| Sopros 617           | CaviWipes™ HP         | No information              |
| Sopros 617           | CaviWipes™ Bleach     | No information              |
| Sopros 717           | CaviCide™             | No information              |
| Sopros 717           | CaviWipes™            | No information              |
| Sopros 717           | Envirocide™           | No information              |
| Sopros 717           | CaviCide1™            | No information              |
| Sopros 717           | CaviWipes1™           | No information              |
| Sopros 717           | CaviWipes™ 2.0        | No information              |
| Sopros 717           | CaviWipes™ HP         | No information              |
| Sopros 717           | CaviWipes™ Bleach     | No information              |
| Sopix Series         | CaviCide™             | Compatible                  |

| <b>Equipment</b> | <b>Metrex Product</b> | <b>Compatibility Status</b> |
|------------------|-----------------------|-----------------------------|
| Sopix Series     | CaviWipes™            | Compatible                  |
| Sopix Series     | Envirocide™           | Compatible                  |
| Sopix Series     | CaviCide1™            | No information              |
| Sopix Series     | CaviWipes1™           | No information              |
| Sopix Series     | CaviWipes™ 2.0        | No information              |
| Sopix Series     | CaviWipes™ HP         | No information              |
| Sopix Series     | CaviWipes™ Bleach     | No information              |
| PSPIX            | CaviCide™             | No information              |
| PSPIX            | CaviWipes™            | No information              |
| PSPIX            | Envirocide™           | No information              |
| PSPIX            | CaviCide1™            | No information              |
| PSPIX            | CaviWipes1™           | No information              |
| PSPIX            | CaviWipes™ 2.0        | No information              |
| PSPIX            | CaviWipes™ HP         | No information              |
| PSPIX            | CaviWipes™ Bleach     | No information              |
| X-Mind AC        | CaviCide™             | No information              |
| X-Mind AC        | CaviWipes™            | No information              |
| X-Mind AC        | Envirocide™           | No information              |
| X-Mind AC        | CaviCide1™            | No information              |
| X-Mind AC        | CaviWipes1™           | No information              |
| X-Mind AC        | CaviWipes™ 2.0        | No information              |
| X-Mind AC        | CaviWipes™ HP         | No information              |
| X-Mind AC        | CaviWipes™ Bleach     | No information              |
| X-Mind Prime     | CaviCide™             | No information              |
| X-Mind Prime     | CaviWipes™            | No information              |
| X-Mind Prime     | Envirocide™           | No information              |

| <b>Equipment</b> | <b>Metrex Product</b> | <b>Compatibility Status</b> |
|------------------|-----------------------|-----------------------------|
| X-Mind Prime     | CaviCide1™            | No information              |
| X-Mind Prime     | CaviWipes1™           | No information              |
| X-Mind Prime     | CaviWipes™ 2.0        | No information              |
| X-Mind Prime     | CaviWipes™ HP         | No information              |
| X-Mind Prime     | CaviWipes™ Bleach     | No information              |
| X-Mind Prime 3D  | CaviCide™             | No information              |
| X-Mind Prime 3D  | CaviWipes™            | No information              |
| X-Mind Prime 3D  | Envirocide™           | No information              |
| X-Mind Prime 3D  | CaviCide1™            | No information              |
| X-Mind Prime 3D  | CaviWipes1™           | No information              |
| X-Mind Prime 3D  | CaviWipes™ 2.0        | No information              |
| X-Mind Prime 3D  | CaviWipes™ HP         | No information              |
| X-Mind Prime 3D  | CaviWipes™ Bleach     | No information              |
| X-Trium          | CaviCide™             | No information              |
| X-Trium          | CaviWipes™            | No information              |
| X-Trium          | Envirocide™           | No information              |
| X-Trium          | CaviCide1™            | No information              |
| X-Trium          | CaviWipes1™           | No information              |
| X-Trium          | CaviWipes™ 2.0        | No information              |
| X-Trium          | CaviWipes™ HP         | No information              |
| X-Trium          | CaviWipes™ Bleach     | No information              |

| <b>Draeger</b>             |                       |                             |
|----------------------------|-----------------------|-----------------------------|
| <b>Equipment</b>           | <b>Metrex Product</b> | <b>Compatibility Status</b> |
| Air-Shield Isolette C400QT | CaviCide™             | Compatible                  |
| Air-Shield Isolette C400QT | CaviWipes™            | Compatible                  |

| <b>Equipment</b>           | <b>Metrex Product</b> | <b>Compatibility Status</b> |
|----------------------------|-----------------------|-----------------------------|
| Air-Shield Isolette C400QT | Envirocide™           | Compatible                  |
| Air-Shield Isolette C400QT | CaviCide1™            | Compatible                  |
| Air-Shield Isolette C400QT | CaviWipes1™           | Compatible                  |
| Air-Shield Isolette C400QT | CaviWipes™ 2.0        | Compatible                  |
| Air-Shield Isolette C400QT | CaviWipes™ HP         | Compatible                  |
| Air-Shield Isolette C400QT | CaviWipes™ Bleach     | Compatible                  |
| Air-Shield Isolette 450QT  | CaviCide™             | Compatible                  |
| Air-Shield Isolette 450QT  | CaviWipes™            | Compatible                  |
| Air-Shield Isolette 450QT  | Envirocide™           | Compatible                  |
| Air-Shield Isolette 450QT  | CaviCide1™            | Compatible                  |
| Air-Shield Isolette 450QT  | CaviWipes1™           | Compatible                  |
| Air-Shield Isolette 450QT  | CaviWipes™ 2.0        | Compatible                  |
| Air-Shield Isolette 450QT  | CaviWipes™ HP         | Compatible                  |
| Air-Shield Isolette 450QT  | CaviWipes™ Bleach     | Compatible                  |
| BabyLeo TN500              | CaviCide™             | No information              |
| BabyLeo TN500              | CaviWipes™            | No information              |
| BabyLeo TN500              | Envirocide™           | No information              |
| BabyLeo TN500              | CaviCide1™            | No information              |
| BabyLeo TN500              | CaviWipes1™           | No information              |
| BabyLeo TN500              | CaviWipes™ 2.0        | No information              |
| BabyLeo TN500              | CaviWipes™ HP         | No information              |
| BabyLeo TN500              | CaviWipes™ Bleach     | No information              |
| Babylog 8000 Plus          | CaviCide™             | Compatible                  |
| Babylog 8000 Plus          | CaviWipes™            | Compatible                  |
| Babylog 8000 Plus          | Envirocide™           | Compatible                  |
| Babylog 8000 Plus          | CaviCide1™            | Compatible                  |

| <b>Equipment</b>                      | <b>Metrex Product</b> | <b>Compatibility Status</b> |
|---------------------------------------|-----------------------|-----------------------------|
| Babylog 8000 Plus                     | CaviWipes1™           | Compatible                  |
| Babylog 8000 Plus                     | CaviWipes™ 2.0        | Compatible                  |
| Babylog 8000 Plus                     | CaviWipes™ HP         | No information              |
| Babylog 8000 Plus                     | CaviWipes™ Bleach     | No information              |
| Babylog VN 500                        | CaviCide™             | No information              |
| Babylog VN 500                        | CaviWipes™            | No information              |
| Babylog VN 500                        | Envirocide™           | No information              |
| Babylog VN 500                        | CaviCide1™            | No information              |
| Babylog VN 500                        | CaviWipes1™           | No information              |
| Babylog VN 500                        | CaviWipes™ 2.0        | No information              |
| Babylog VN 500                        | CaviWipes™ HP         | No information              |
| Babylog VN 500                        | CaviWipes™ Bleach     | No information              |
| D-Vapor                               | CaviCide™             | No information              |
| D-Vapor                               | CaviWipes™            | No information              |
| D-Vapor                               | Envirocide™           | No information              |
| D-Vapor                               | CaviCide1™            | No information              |
| Fabius MRI all respiratory components | CaviWipes1™           | Compatible                  |
| Fabius MRI all respiratory components | CaviWipes™ 2.0        | Compatible                  |
| Fabius MRI all respiratory components | CaviWipes™ HP         | Not recommended for use     |
| Fabius MRI all respiratory components | CaviWipes™ Bleach     | Compatible                  |
| Fabius MRI Spirolog Flow Sensors      | CaviCide™             | No information              |
| Fabius MRI Spirolog Flow Sensors      | CaviWipes™            | No information              |

| <b>Equipment</b>                 | <b>Metrex Product</b> | <b>Compatibility Status</b> |
|----------------------------------|-----------------------|-----------------------------|
| Fabius MRI Spirolog Flow Sensors | Envirocide™           | No information              |
| Fabius MRI Spirolog Flow Sensors | CaviCide1™            | No information              |
| Fabius MRI Spirolog Flow Sensors | CaviWipes1™           | No information              |
| Fabius MRI Spirolog Flow Sensors | CaviWipes™ 2.0        | No information              |
| Fabius MRI Spirolog Flow Sensors | CaviWipes™ HP         | No information              |
| Fabius MRI Spirolog Flow Sensors | CaviWipes™ Bleach     | No information              |
| Fabius MRI surface components    | CaviCide™             | No information              |
| Fabius MRI surface components    | CaviWipes™            | No information              |
| Fabius MRI surface components    | Envirocide™           | No information              |
| Fabius MRI surface components    | CaviCide1™            | No information              |
| Fabius MRI surface components    | CaviWipes1™           | No information              |
| Fabius MRI surface components    | CaviWipes™ 2.0        | No information              |
| Fabius MRI surface components    | CaviWipes™ HP         | No information              |
| Fabius MRI surface components    | CaviWipes™ Bleach     | No information              |
| Fabius Plus XL Workstation       | CaviCide™             | No information              |
| Fabius Plus XL Workstation       | CaviWipes™            | No information              |
| Fabius Plus XL Workstation       | Envirocide™           | No information              |



| <b>Equipment</b>                     | <b>Metrex Product</b> | <b>Compatibility Status</b> |
|--------------------------------------|-----------------------|-----------------------------|
| Fabius Plus XL Workstation           | CaviCide1™            | No information              |
| Fabius Plus XL Workstation           | CaviWipes1™           | No information              |
| Fabius Plus XL Workstation           | CaviWipes™ 2.0        | No information              |
| Fabius Plus XL Workstation           | CaviWipes™ HP         | No information              |
| Fabius Plus XL Workstation           | CaviWipes™ Bleach     | No information              |
| Fabius Trio non-critical components  | CaviCide™             | No information              |
| Fabius Trio non-critical components  | CaviWipes™            | No information              |
| Fabius Trio non-critical components  | Envirocide™           | No information              |
| Fabius Trio non-critical components  | CaviCide1™            | No information              |
| Fabius Trio non-critical components  | CaviWipes1™           | No information              |
| Fabius Trio non-critical components  | CaviWipes™ 2.0        | No information              |
| Fabius Trio non-critical components  | CaviWipes™ HP         | No information              |
| Fabius Trio non-critical components  | CaviWipes™ Bleach     | No information              |
| Fabius Trio semi-critical components | CaviCide™             | No information              |
| Fabius Trio semi-critical components | CaviWipes™            | No information              |
| Fabius Trio semi-critical components | Envirocide™           | No information              |
| Fabius Trio semi-critical components | CaviCide1™            | No information              |
| D-Vapor                              | CaviWipes1™           | No information              |

| <b>Equipment</b>    | <b>Metrex Product</b> | <b>Compatibility Status</b> |
|---------------------|-----------------------|-----------------------------|
| D-Vapor             | CaviWipes™ 2.0        | No information              |
| D-Vapor             | CaviWipes™ HP         | No information              |
| D-Vapor             | CaviWipes™ Bleach     | No information              |
| D-Vapor 3000        | CaviCide™             | No information              |
| D-Vapor 3000        | CaviWipes™            | No information              |
| D-Vapor 3000        | Envirocide™           | No information              |
| D-Vapor 3000        | CaviCide1™            | No information              |
| D-Vapor 3000        | CaviWipes1™           | No information              |
| D-Vapor 3000        | CaviWipes™ 2.0        | No information              |
| D-Vapor 3000        | CaviWipes™ HP         | No information              |
| D-Vapor 3000        | CaviWipes™ Bleach     | No information              |
| Evita Infinity V500 | CaviCide™             | No information              |
| Evita Infinity V500 | CaviWipes™            | No information              |
| Evita Infinity V500 | Envirocide™           | No information              |
| Evita Infinity V500 | CaviCide1™            | No information              |
| Evita Infinity V500 | CaviWipes1™           | No information              |
| Evita Infinity V500 | CaviWipes™ 2.0        | No information              |
| Evita Infinity V500 | CaviWipes™ HP         | No information              |
| Evita Infinity V500 | CaviWipes™ Bleach     | No information              |
| Evita XL            | CaviCide™             | Compatible                  |
| Evita XL            | CaviWipes™            | Compatible                  |
| Evita XL            | Envirocide™           | Compatible                  |
| Evita XL            | CaviCide1™            | Compatible                  |
| Evita XL            | CaviWipes1™           | Compatible                  |
| Evita XL            | CaviWipes™ 2.0        | Compatible                  |
| Evita XL            | CaviWipes™ HP         | No information              |

| <b>Equipment</b>                      | <b>Metrex Product</b> | <b>Compatibility Status</b> |
|---------------------------------------|-----------------------|-----------------------------|
| Evita XL                              | CaviWipes™ Bleach     | No information              |
| Fabius GS Premium                     | CaviCide™             | No information              |
| Fabius GS Premium                     | CaviWipes™            | No information              |
| Fabius GS Premium                     | Envirocide™           | No information              |
| Fabius GS Premium                     | CaviCide1™            | No information              |
| Fabius GS Premium                     | CaviWipes1™           | No information              |
| Fabius GS Premium                     | CaviWipes™ 2.0        | No information              |
| Fabius GS Premium                     | CaviWipes™ HP         | No information              |
| Fabius GS Premium                     | CaviWipes™ Bleach     | No information              |
| Fabius MRI all respiratory components | CaviCide™             | Compatible                  |
| Fabius MRI all respiratory components | CaviWipes™            | Compatible                  |
| Fabius MRI all respiratory components | Envirocide™           | Compatible                  |
| Fabius MRI all respiratory components | CaviCide1™            | Compatible                  |
| Savina 300                            | CaviWipes1™           | No information              |
| Savina 300                            | CaviWipes™ 2.0        | No information              |
| Savina 300                            | CaviWipes™ HP         | No information              |
| Savina 300                            | CaviWipes™ Bleach     | No information              |
| VN Series                             | CaviCide™             | No information              |
| VN Series                             | CaviWipes™            | No information              |
| VN Series                             | Envirocide™           | No information              |
| VN Series                             | CaviCide1™            | No information              |
| VN Series                             | CaviWipes1™           | No information              |
| VN Series                             | CaviWipes™ 2.0        | No information              |

| <b>Equipment</b> | <b>Metrex Product</b> | <b>Compatibility Status</b> |
|------------------|-----------------------|-----------------------------|
| VN Series        | CaviWipes™ HP         | No information              |
| VN Series        | CaviWipes™ Bleach     | No information              |
| X-am 3500        | CaviCide™             | No information              |
| X-am 3500        | CaviWipes™            | No information              |
| X-am 3500        | Envirocide™           | No information              |
| X-am 3500        | CaviCide1™            | No information              |
| X-am 3500        | CaviWipes1™           | No information              |
| X-am 3500        | CaviWipes™ 2.0        | No information              |
| X-am 3500        | CaviWipes™ HP         | No information              |
| X-am 3500        | CaviWipes™ Bleach     | No information              |
| X-am 5000        | CaviCide™             | No information              |
| X-am 5000        | CaviWipes™            | No information              |
| X-am 5000        | Envirocide™           | No information              |
| X-am 5000        | CaviCide1™            | No information              |
| X-am 5000        | CaviWipes1™           | No information              |
| X-am 5000        | CaviWipes™ 2.0        | No information              |
| X-am 5000        | CaviWipes™ HP         | No information              |
| X-am 5000        | CaviWipes™ Bleach     | No information              |
| X-Dock 6300      | CaviCide™             | No information              |
| X-Dock 6300      | CaviWipes™            | No information              |
| X-Dock 6300      | Envirocide™           | No information              |
| X-Dock 6300      | CaviCide1™            | No information              |
| X-Dock 6300      | CaviWipes1™           | No information              |
| X-Dock 6300      | CaviWipes™ 2.0        | No information              |
| X-Dock 6300      | CaviWipes™ HP         | No information              |
| X-Dock 6300      | CaviWipes™ Bleach     | No information              |

| <b>Equipment</b>                     | <b>Metrex Product</b> | <b>Compatibility Status</b> |
|--------------------------------------|-----------------------|-----------------------------|
| Fabius Trio semi-critical components | CaviWipes1™           | No information              |
| Fabius Trio semi-critical components | CaviWipes™ 2.0        | No information              |
| Fabius Trio semi-critical components | CaviWipes™ HP         | No information              |
| Fabius Trio semi-critical components | CaviWipes™ Bleach     | No information              |
| Oxylog 3000 Plus                     | CaviCide™             | No information              |
| Oxylog 3000 Plus                     | CaviWipes™            | No information              |
| Oxylog 3000 Plus                     | Envirocide™           | No information              |
| Oxylog 3000 Plus                     | CaviCide1™            | No information              |
| Oxylog 3000 Plus                     | CaviWipes1™           | No information              |
| Oxylog 3000 Plus                     | CaviWipes™ 2.0        | No information              |
| Oxylog 3000 Plus                     | CaviWipes™ HP         | No information              |
| Oxylog 3000 Plus                     | CaviWipes™ Bleach     | No information              |
| Perseus A500 non-critical components | CaviCide™             | No information              |
| Perseus A500 non-critical components | CaviWipes™            | No information              |
| Perseus A500 non-critical components | Envirocide™           | No information              |
| Perseus A500 non-critical components | CaviCide1™            | No information              |
| Perseus A500 non-critical components | CaviWipes1™           | No information              |
| Perseus A500 non-critical components | CaviWipes™ 2.0        | No information              |
| Perseus A500 non-critical components | CaviWipes™ HP         | No information              |

| <b>Equipment</b>                      | <b>Metrex Product</b> | <b>Compatibility Status</b> |
|---------------------------------------|-----------------------|-----------------------------|
| Perseus A500 non-critical components  | CaviWipes™ Bleach     | No information              |
| Perseus A500 semi-critical components | CaviCide™             | No information              |
| Perseus A500 semi-critical components | CaviWipes™            | No information              |
| Perseus A500 semi-critical components | Envirocide™           | No information              |
| Perseus A500 semi-critical components | CaviCide1™            | No information              |
| Perseus A500 semi-critical components | CaviWipes1™           | No information              |
| Perseus A500 semi-critical components | CaviWipes™ 2.0        | No information              |
| Perseus A500 semi-critical components | CaviWipes™ HP         | No information              |
| Perseus A500 semi-critical components | CaviWipes™ Bleach     | No information              |
| Primus                                | CaviCide™             | No information              |
| Primus                                | CaviWipes™            | No information              |
| Primus                                | Envirocide™           | No information              |
| Primus                                | CaviCide1™            | No information              |
| Primus                                | CaviWipes1™           | No information              |
| Primus                                | CaviWipes™ 2.0        | No information              |
| Primus                                | CaviWipes™ HP         | No information              |
| Primus                                | CaviWipes™ Bleach     | No information              |
| Savina 300                            | CaviCide™             | No information              |
| Savina 300                            | CaviWipes™            | No information              |
| Savina 300                            | Envirocide™           | No information              |

| <b>Equipment</b> | <b>Metrex Product</b> | <b>Compatibility Status</b> |
|------------------|-----------------------|-----------------------------|
| Savina 300       | CaviCide1™            | No information              |

| <b>KMI</b>             |                       |                             |
|------------------------|-----------------------|-----------------------------|
| <b>Equipment</b>       | <b>Metrex Product</b> | <b>Compatibility Status</b> |
| Phantom 2550-1 SUCTION | CaviCide™             | No information              |
| Phantom 2550-1 SUCTION | CaviWipes™            | No information              |
| Phantom 2550-1 SUCTION | Envirocide™           | No information              |
| Phantom 2550-1 SUCTION | CaviCide1™            | No information              |
| Phantom 2550-1 SUCTION | CaviWipes1™           | No information              |
| Phantom 2550-1 SUCTION | CaviWipes™ 2.0        | No information              |
| Phantom 2550-1 SUCTION | CaviWipes™ HP         | No information              |
| Phantom 2550-1 SUCTION | CaviWipes™ Bleach     | No information              |
| Phantom 2550-1 ASP     | CaviCide™             | No information              |
| Phantom 2550-1 ASP     | CaviWipes™            | No information              |
| Phantom 2550-1 ASP     | Envirocide™           | No information              |
| Phantom 2550-1 ASP     | CaviCide1™            | No information              |
| Phantom 2550-1 ASP     | CaviWipes1™           | No information              |
| Phantom 2550-1 ASP     | CaviWipes™ 2.0        | No information              |
| Phantom 2550-1 ASP     | CaviWipes™ HP         | No information              |
| Phantom 2550-1 ASP     | CaviWipes™ Bleach     | No information              |

| <b>Roche</b>     |                       |                             |
|------------------|-----------------------|-----------------------------|
| <b>Equipment</b> | <b>Metrex Product</b> | <b>Compatibility Status</b> |
| Cobas Plus       | CaviWipes1™           | Compatible                  |
| Cobas Plus       | CaviWipes™ 2.0        | Compatible                  |
| Cobas Plus       | CaviWipes™ HP         | No information              |

| <b>Equipment</b>    | <b>Metrex Product</b> | <b>Compatibility Status</b> |
|---------------------|-----------------------|-----------------------------|
| Cobas Plus          | CaviWipes™ Bleach     | Compatible                  |
| cobas u 601         | CaviCide™             | No information              |
| cobas u 601         | CaviWipes™            | No information              |
| cobas u 601         | Envirocide™           | No information              |
| cobas u 601         | CaviCide1™            | No information              |
| cobas u 601         | CaviWipes1™           | No information              |
| cobas u 601         | CaviWipes™ 2.0        | No information              |
| cobas u 601         | CaviWipes™ HP         | No information              |
| cobas u 601         | CaviWipes™ Bleach     | No information              |
| Cobas Urisys 1100   | CaviCide™             | No information              |
| Cobas Urisys 1100   | CaviWipes™            | No information              |
| Cobas Urisys 1100   | Envirocide™           | No information              |
| Cobas Urisys 1100   | CaviCide1™            | No information              |
| Cobas Urisys 1100   | CaviWipes1™           | No information              |
| Cobas Urisys 1100   | CaviWipes™ 2.0        | No information              |
| Cobas Urisys 1100   | CaviWipes™ HP         | No information              |
| Cobas Urisys 1100   | CaviWipes™ Bleach     | No information              |
| Model 9180          | CaviCide™             | No information              |
| Model 9180          | CaviWipes™            | No information              |
| Model 9180          | Envirocide™           | No information              |
| Model 9180          | CaviCide1™            | No information              |
| Model 9180          | CaviWipes1™           | No information              |
| Model 9180          | CaviWipes™ 2.0        | No information              |
| Model 9180          | CaviWipes™ HP         | No information              |
| Model 9180          | CaviWipes™ Bleach     | No information              |
| Multiplate Analyzer | CaviCide™             | No information              |
| Multiplate Analyzer | CaviWipes™            | No information              |



| <b>Equipment</b>                | <b>Metrex Product</b> | <b>Compatibility Status</b> |
|---------------------------------|-----------------------|-----------------------------|
| Multiplate Analyzer             | Envirocide™           | No information              |
| Multiplate Analyzer             | CaviCide1™            | No information              |
| Multiplate Analyzer             | CaviWipes1™           | No information              |
| Multiplate Analyzer             | CaviWipes™ 2.0        | No information              |
| Multiplate Analyzer             | CaviWipes™ HP         | No information              |
| Multiplate Analyzer             | CaviWipes™ Bleach     | No information              |
| CoaguChek Vantus                | CaviCide™             | No information              |
| CoaguChek Vantus                | CaviWipes™            | No information              |
| CoaguChek Vantus                | Envirocide™           | No information              |
| CoaguChek Vantus                | CaviCide1™            | No information              |
| CoaguChek Vantus                | CaviWipes1™           | No information              |
| CoaguChek Vantus                | CaviWipes™ 2.0        | No information              |
| CoaguChek Vantus                | CaviWipes™ HP         | No information              |
| CoaguChek Vantus                | CaviWipes™ Bleach     | No information              |
| CoaguChek XS Plus/meter housing | CaviCide™             | No information              |
| CoaguChek XS Plus/meter housing | CaviWipes™            | No information              |
| CoaguChek XS Plus/meter housing | Envirocide™           | No information              |
| CoaguChek XS Plus/meter housing | CaviCide1™            | No information              |
| CoaguChek XS Plus/meter housing | CaviWipes1™           | No information              |
| CoaguChek XS Plus/meter housing | CaviWipes™ 2.0        | No information              |
| CoaguChek XS Plus/meter housing | CaviWipes™ HP         | No information              |
| CoaguChek XS Plus/meter housing | CaviWipes™ Bleach     | Compatible                  |

| <b>Equipment</b>                      | <b>Metrex Product</b> | <b>Compatibility Status</b> |
|---------------------------------------|-----------------------|-----------------------------|
| CoaguChek XS Pro                      | CaviCide™             | No information              |
| CoaguChek XS Pro                      | CaviWipes™            | No information              |
| CoaguChek XS Pro                      | Envirocide™           | No information              |
| CoaguChek XS Pro                      | CaviCide1™            | No information              |
| CoaguChek XS Pro                      | CaviWipes1™           | No information              |
| CoaguChek XS Pro                      | CaviWipes™ 2.0        | No information              |
| CoaguChek XS Pro                      | CaviWipes™ HP         | No information              |
| CoaguChek XS Pro                      | CaviWipes™ Bleach     | No information              |
| Coagu-Chek XS System                  | CaviCide™             | No information              |
| Coagu-Chek XS System                  | CaviWipes™            | No information              |
| Coagu-Chek XS System                  | Envirocide™           | No information              |
| Coagu-Chek XS System                  | CaviCide1™            | No information              |
| Coagu-Chek XS System                  | CaviWipes1™           | No information              |
| Coagu-Chek XS System                  | CaviWipes™ 2.0        | No information              |
| Coagu-Chek XS System                  | CaviWipes™ HP         | No information              |
| Coagu-Chek XS System                  | CaviWipes™ Bleach     | No information              |
| Coagu-Chek XS System-TEST STRIP GUIDE | CaviCide™             | No information              |
| Coagu-Chek XS System-TEST STRIP GUIDE | CaviWipes™            | No information              |
| Coagu-Chek XS System-TEST STRIP GUIDE | Envirocide™           | No information              |
| Coagu-Chek XS System-TEST STRIP GUIDE | CaviCide1™            | No information              |
| Coagu-Chek XS System-TEST STRIP GUIDE | CaviWipes1™           | No information              |
| Coagu-Chek XS System-TEST STRIP GUIDE | CaviWipes™ 2.0        | No information              |

| <b>Equipment</b>                      | <b>Metrex Product</b> | <b>Compatibility Status</b> |
|---------------------------------------|-----------------------|-----------------------------|
| Coagu-Chek XS System-TEST STRIP GUIDE | CaviWipes™ HP         | No information              |
| Coagu-Chek XS System-TEST STRIP GUIDE | CaviWipes™ Bleach     | Compatible                  |
| Accutrend Plus                        | CaviWipes1™           | No information              |
| Accutrend Plus                        | CaviWipes™ 2.0        | No information              |
| Accutrend Plus                        | CaviWipes™ HP         | No information              |
| Accutrend Plus                        | CaviWipes™ Bleach     | No information              |
| CoaguChek InRange                     | CaviCide™             | No information              |
| CoaguChek InRange                     | CaviWipes™            | No information              |
| CoaguChek InRange                     | Envirocide™           | No information              |
| CoaguChek InRange                     | CaviCide1™            | No information              |
| CoaguChek InRange                     | CaviWipes1™           | No information              |
| CoaguChek InRange                     | CaviWipes™ 2.0        | No information              |
| CoaguChek InRange                     | CaviWipes™ HP         | No information              |
| CoaguChek InRange                     | CaviWipes™ Bleach     | No information              |
| CoaguChek Pro II                      | CaviCide™             | Compatible                  |
| CoaguChek Pro II                      | CaviWipes™            | Compatible                  |
| CoaguChek Pro II                      | Envirocide™           | Compatible                  |
| CoaguChek Pro II                      | CaviCide1™            | Compatible                  |
| CoaguChek Pro II                      | CaviWipes1™           | Compatible                  |
| CoaguChek Pro II                      | CaviWipes™ 2.0        | Compatible                  |
| CoaguChek Pro II                      | CaviWipes™ HP         | No information              |
| CoaguChek Pro II                      | CaviWipes™ Bleach     | Compatible                  |
| Cobas 4800 System                     | CaviCide™             | No information              |
| Cobas 4800 System                     | CaviWipes™            | No information              |
| Cobas 4800 System                     | Envirocide™           | No information              |
| Cobas 4800 System                     | CaviCide1™            | No information              |

| <b>Equipment</b>               | <b>Metrex Product</b> | <b>Compatibility Status</b> |
|--------------------------------|-----------------------|-----------------------------|
| Cobas 4800 System              | CaviWipes1™           | No information              |
| Cobas 4800 System              | CaviWipes™ 2.0        | No information              |
| Cobas 4800 System              | CaviWipes™ HP         | No information              |
| Cobas 4800 System              | CaviWipes™ Bleach     | No information              |
| Cobas Miras systems            | CaviCide™             | No information              |
| Cobas Miras systems            | CaviWipes™            | No information              |
| Cobas Miras systems            | Envirocide™           | No information              |
| Cobas Miras systems            | CaviCide1™            | No information              |
| Cobas Miras systems            | CaviWipes1™           | No information              |
| Cobas Miras systems            | CaviWipes™ 2.0        | No information              |
| Cobas Miras systems            | CaviWipes™ HP         | No information              |
| Cobas Miras systems            | CaviWipes™ Bleach     | Compatible                  |
| Cobas Plus                     | CaviCide™             | Compatible                  |
| Cobas Plus                     | CaviWipes™            | Compatible                  |
| Cobas Plus                     | Envirocide™           | Compatible                  |
| Cobas Plus                     | CaviCide1™            | Compatible                  |
| Accu-Check Inform II           | CaviCide™             | No information              |
| Accu-Check Inform II           | CaviWipes™            | No information              |
| Accu-Check Inform II           | Envirocide™           | No information              |
| Accu-Check Inform II           | CaviCide1™            | No information              |
| Accu-Check Inform II           | CaviWipes1™           | No information              |
| Accu-Check Inform II           | CaviWipes™ 2.0        | No information              |
| Accu-Check Inform II           | CaviWipes™ HP         | No information              |
| Accu-Check Inform II           | CaviWipes™ Bleach     | No information              |
| Accu-check Inform II Base Unit | CaviCide™             | No information              |
| Accu-check Inform II Base Unit | CaviWipes™            | No information              |

| <b>Equipment</b>               | <b>Metrex Product</b> | <b>Compatibility Status</b> |
|--------------------------------|-----------------------|-----------------------------|
| Accu-check Inform II Base Unit | Envirocide™           | No information              |
| Accu-check Inform II Base Unit | CaviCide1™            | No information              |
| Accu-check Inform II Base Unit | CaviWipes1™           | No information              |
| Accu-check Inform II Base Unit | CaviWipes™ 2.0        | No information              |
| Accu-check Inform II Base Unit | CaviWipes™ HP         | No information              |
| Accu-check Inform II Base Unit | CaviWipes™ Bleach     | No information              |
| Accutrend Plus                 | CaviCide™             | No information              |
| Accutrend Plus                 | CaviWipes™            | No information              |
| Accutrend Plus                 | Envirocide™           | No information              |
| Accutrend Plus                 | CaviCide1™            | No information              |