

# Material Compatibility Status

| A-dec                                              |                   |                         |
|----------------------------------------------------|-------------------|-------------------------|
| Equipment                                          | Metrex Product    | Compatibility Status    |
| Preference Collection Model 5562 Laminate & quartz | CaviWipes1™       | No information          |
| Preference Collection Model 5562 Laminate & quartz | CaviWipes™ 2.0    | No information          |
| Preference Collection Model 5562 Laminate & quartz | CaviWipes™ HP     | No information          |
| Preference Collection Model 5562 Laminate & quartz | CaviWipes™ Bleach | Not recommended for use |
| Preference Collection Model 5580 Laminate & quartz | CaviCide™         | No information          |
| Preference Collection Model 5580 Laminate & quartz | CaviWipes™        | No information          |
| Preference Collection Model 5580 Laminate & quartz | Envirocide™       | No information          |
| Preference Collection Model 5580 Laminate & quartz | CaviCide1™        | No information          |
| Preference Collection Model 5580 Laminate & quartz | CaviWipes1™       | No information          |
| Preference Collection Model 5580 Laminate & quartz | CaviWipes™ 2.0    | No information          |
| Preference Collection Model 5580 Laminate & quartz | CaviWipes™ HP     | No information          |
| Preference Collection Model 5580 Laminate & quartz | CaviWipes™ Bleach | Not recommended for use |
| Inspire Collection: Infills                        | CaviCide™         | Compatible              |
| Inspire Collection: Infills                        | CaviWipes™        | Compatible              |

| <b>Equipment</b>                                   | <b>Metrex Product</b> | <b>Compatibility Status</b> |
|----------------------------------------------------|-----------------------|-----------------------------|
| Inspire Collection: Infills                        | Envirocide™           | Compatible                  |
| Inspire Collection: Infills                        | CaviCide1™            | Compatible                  |
| Inspire Collection: Infills                        | CaviWipes1™           | Compatible                  |
| Inspire Collection: Infills                        | CaviWipes™ 2.0        | Compatible                  |
| Inspire Collection: Infills                        | CaviWipes™ HP         | No information              |
| Inspire Collection: Infills                        | CaviWipes™ Bleach     | Not recommended for use     |
| Preference Collection Model 5562 Laminate & quartz | CaviCide™             | No information              |
| Preference Collection Model 5562 Laminate & quartz | CaviWipes™            | No information              |
| Preference Collection Model 5562 Laminate & quartz | Envirocide™           | No information              |
| Preference Collection Model 5562 Laminate & quartz | CaviCide1™            | No information              |
| Inspire Collection: Laminate surfaces              | CaviCide™             | Compatible                  |
| Inspire Collection: Laminate surfaces              | CaviWipes™            | Compatible                  |
| Inspire Collection: Laminate surfaces              | Envirocide™           | Compatible                  |
| Inspire Collection: Laminate surfaces              | CaviCide1™            | Compatible                  |
| Inspire Collection: Laminate surfaces              | CaviWipes1™           | Compatible                  |
| Inspire Collection: Laminate surfaces              | CaviWipes™ 2.0        | Compatible                  |
| Inspire Collection: Laminate surfaces              | CaviWipes™ HP         | No information              |
| Inspire Collection: Laminate surfaces              | CaviWipes™ Bleach     | Not recommended for use     |

| <b>Equipment</b>                    | <b>Metrex Product</b> | <b>Compatibility Status</b> |
|-------------------------------------|-----------------------|-----------------------------|
| Inspire Collection: Quartz surfaces | CaviCide™             | Compatible                  |
| Inspire Collection: Quartz surfaces | CaviWipes™            | Compatible                  |
| Inspire Collection: Quartz surfaces | Envirocide™           | Compatible                  |
| Inspire Collection: Quartz surfaces | CaviCide1™            | Compatible                  |
| Inspire Collection: Quartz surfaces | CaviWipes1™           | Compatible                  |
| Inspire Collection: Quartz surfaces | CaviWipes™ 2.0        | Compatible                  |
| Inspire Collection: Quartz surfaces | CaviWipes™ HP         | No information              |
| Inspire Collection: Quartz surfaces | CaviWipes™ Bleach     | Not recommended for use     |
| Inspire Collection: Solid surfaces  | CaviCide™             | Compatible                  |
| Inspire Collection: Solid surfaces  | CaviWipes™            | Compatible                  |
| Inspire Collection: Solid surfaces  | Envirocide™           | Compatible                  |
| Inspire Collection: Solid surfaces  | CaviCide1™            | Compatible                  |
| Inspire Collection: Solid surfaces  | CaviWipes1™           | Compatible                  |
| Inspire Collection: Solid surfaces  | CaviWipes™ 2.0        | Compatible                  |
| Inspire Collection: Solid surfaces  | CaviWipes™ HP         | No information              |
| Inspire Collection: Solid surfaces  | CaviWipes™ Bleach     | Not recommended for use     |

| <b>Equipment</b> | <b>Metrex Product</b> | <b>Compatibility Status</b> |
|------------------|-----------------------|-----------------------------|
| A-dec 400        | CaviCide™             | Compatible                  |
| A-dec 400        | CaviWipes™            | Compatible                  |
| A-dec 400        | Envirocide™           | Compatible                  |
| A-dec 400        | CaviCide1™            | Compatible                  |
| A-dec 400        | CaviWipes1™           | Compatible                  |
| A-dec 400        | CaviWipes™ 2.0        | Compatible                  |
| A-dec 400        | CaviWipes™ HP         | No information              |
| A-dec 400        | CaviWipes™ Bleach     | Not recommended for use     |
| A-dec 500        | CaviCide™             | Compatible                  |
| A-dec 500        | CaviWipes™            | Compatible                  |
| A-dec 500        | Envirocide™           | Compatible                  |
| A-dec 500        | CaviCide1™            | Compatible                  |
| A-dec 500        | CaviWipes1™           | Compatible                  |
| A-dec 500        | CaviWipes™ 2.0        | Compatible                  |
| A-dec 500        | CaviWipes™ HP         | No information              |
| A-dec 500        | CaviWipes™ Bleach     | Not recommended for use     |
| A-dec 511b       | CaviCide™             | Compatible                  |
| A-dec 511b       | CaviWipes™            | Compatible                  |
| A-dec 511b       | Envirocide™           | Compatible                  |
| A-dec 511b       | CaviCide1™            | Compatible                  |
| A-dec 511b       | CaviWipes1™           | Compatible                  |
| A-dec 511b       | CaviWipes™ 2.0        | Compatible                  |
| A-dec 511b       | CaviWipes™ HP         | Not recommended for use     |
| A-dec 511b       | CaviWipes™ Bleach     | Not recommended for use     |
| A-dec 300        | CaviCide™             | Compatible                  |
| A-dec 300        | CaviWipes™            | Compatible                  |
| A-dec 300        | Envirocide™           | Compatible                  |

| <b>Equipment</b>                                                          | <b>Metrex Product</b> | <b>Compatibility Status</b> |
|---------------------------------------------------------------------------|-----------------------|-----------------------------|
| A-dec 300                                                                 | CaviCide1™            | Compatible                  |
| Model 200                                                                 | CaviCide™             | No information              |
| Model 200                                                                 | CaviWipes™            | No information              |
| Model 200                                                                 | Envirocide™           | No information              |
| Model 200                                                                 | CaviCide1™            | No information              |
| Model 200                                                                 | CaviWipes1™           | No information              |
| Model 200                                                                 | CaviWipes™ 2.0        | No information              |
| Model 200                                                                 | CaviWipes™ HP         | No information              |
| Model 200                                                                 | CaviWipes™ Bleach     | No information              |
| A-dec 300                                                                 | CaviWipes1™           | Compatible                  |
| A-dec 300                                                                 | CaviWipes™ 2.0        | Compatible                  |
| A-dec 300                                                                 | CaviWipes™ HP         | No information              |
| A-dec 300                                                                 | CaviWipes™ Bleach     | Not recommended for use     |
| Preference Collection Model 5562 Treatment console insert / solid surface | CaviCide™             | No information              |
| Preference Collection Model 5562 Treatment console insert / solid surface | CaviWipes™            | No information              |
| Preference Collection Model 5562 Treatment console insert / solid surface | Envirocide™           | No information              |
| Preference Collection Model 5562 Treatment console insert / solid surface | CaviCide1™            | No information              |
| Preference Collection Model 5562 Treatment console insert / solid surface | CaviWipes1™           | No information              |
| Preference Collection Model 5562 Treatment console insert / solid surface | CaviWipes™ 2.0        | No information              |

| <b>Equipment</b>                                                          | <b>Metrex Product</b> | <b>Compatibility Status</b> |
|---------------------------------------------------------------------------|-----------------------|-----------------------------|
| Preference Collection Model 5562 Treatment console insert / solid surface | CaviWipes™ HP         | No information              |
| Preference Collection Model 5562 Treatment console insert / solid surface | CaviWipes™ Bleach     | Not recommended for use     |
| Preference Collection Model 5580 Treatment console / solid surface        | CaviCide™             | No information              |
| Preference Collection Model 5580 Treatment console / solid surface        | CaviWipes™            | No information              |
| Preference Collection Model 5580 Treatment console / solid surface        | Envirocide™           | No information              |
| Preference Collection Model 5580 Treatment console / solid surface        | CaviCide1™            | No information              |
| Preference Collection Model 5580 Treatment console / solid surface        | CaviWipes1™           | No information              |
| Preference Collection Model 5580 Treatment console / solid surface        | CaviWipes™ 2.0        | No information              |
| Preference Collection Model 5580 Treatment console / solid surface        | CaviWipes™ HP         | No information              |
| Preference Collection Model 5580 Treatment console / solid surface        | CaviWipes™ Bleach     | Not recommended for use     |
| Inspire 391 Treatment Console                                             | CaviCide™             | Compatible                  |
| Inspire 391 Treatment Console                                             | CaviWipes™            | Compatible                  |

| <b>Equipment</b>                | <b>Metrex Product</b> | <b>Compatibility Status</b> |
|---------------------------------|-----------------------|-----------------------------|
| Inspire 391 Treatment Console   | Envirocide™           | Compatible                  |
| Inspire 391 Treatment Console   | CaviCide1™            | Compatible                  |
| Inspire 391 Treatment Console   | CaviWipes1™           | Compatible                  |
| Inspire 391 Treatment Console   | CaviWipes™ 2.0        | Compatible                  |
| Inspire 391 Treatment Console   | CaviWipes™ HP         | No information              |
| Inspire 391 Treatment Console   | CaviWipes™ Bleach     | Not recommended for use     |
| Inspire 591 A Treatment Console | CaviCide™             | Compatible                  |
| Inspire 591 A Treatment Console | CaviWipes™            | Compatible                  |
| Inspire 591 A Treatment Console | Envirocide™           | Compatible                  |
| Inspire 591 A Treatment Console | CaviCide1™            | Compatible                  |
| Inspire 591 A Treatment Console | CaviWipes1™           | Compatible                  |
| Inspire 591 A Treatment Console | CaviWipes™ 2.0        | Compatible                  |
| Inspire 591 A Treatment Console | CaviWipes™ HP         | No information              |
| Inspire 591 A Treatment Console | CaviWipes™ Bleach     | Not recommended for use     |
| Inspire 592 A Central Console   | CaviCide™             | Compatible                  |
| Inspire 592 A Central Console   | CaviWipes™            | Compatible                  |

| <b>Equipment</b>              | <b>Metrex Product</b> | <b>Compatibility Status</b> |
|-------------------------------|-----------------------|-----------------------------|
| Inspire 592 A Central Console | Envirocide™           | Compatible                  |
| Inspire 592 A Central Console | CaviCide1™            | Compatible                  |
| Inspire 592 A Central Console | CaviWipes1™           | Compatible                  |
| Inspire 592 A Central Console | CaviWipes™ 2.0        | Compatible                  |
| Inspire 592 A Central Console | CaviWipes™ HP         | No information              |
| Inspire 592 A Central Console | CaviWipes™ Bleach     | Not recommended for use     |
| Inspire 593 a Side Console    | CaviCide™             | Compatible                  |
| Inspire 593 a Side Console    | CaviWipes™            | Compatible                  |
| Inspire 593 a Side Console    | Envirocide™           | Compatible                  |
| Inspire 593 a Side Console    | CaviCide1™            | Compatible                  |
| A-dec 311                     | CaviCide™             | Compatible                  |
| A-dec 311                     | CaviWipes™            | Compatible                  |
| A-dec 311                     | Envirocide™           | Compatible                  |
| A-dec 311                     | CaviCide1™            | Compatible                  |
| A-dec 311                     | CaviWipes1™           | Compatible                  |
| A-dec 311                     | CaviWipes™ 2.0        | Compatible                  |
| A-dec 311                     | CaviWipes™ HP         | Not recommended for use     |
| A-dec 311                     | CaviWipes™ Bleach     | Not recommended for use     |
| A-dec 411                     | CaviCide™             | Compatible                  |
| A-dec 411                     | CaviWipes™            | Compatible                  |
| A-dec 411                     | Envirocide™           | Compatible                  |
| A-dec 411                     | CaviCide1™            | Compatible                  |
| A-dec 411                     | CaviWipes1™           | Compatible                  |

| <b>Equipment</b>                                       | <b>Metrex Product</b> | <b>Compatibility Status</b> |
|--------------------------------------------------------|-----------------------|-----------------------------|
| A-dec 411                                              | CaviWipes™ 2.0        | Compatible                  |
| A-dec 411                                              | CaviWipes™ HP         | Not recommended for use     |
| A-dec 411                                              | CaviWipes™ Bleach     | Not recommended for use     |
| Inspire 593 a Side Console                             | CaviWipes1™           | Compatible                  |
| Inspire 593 a Side Console                             | CaviWipes™ 2.0        | Compatible                  |
| Inspire 593 a Side Console                             | CaviWipes™ HP         | No information              |
| Inspire 593 a Side Console                             | CaviWipes™ Bleach     | Not recommended for use     |
| Inspire 594 A Sterilization Center                     | CaviCide™             | Compatible                  |
| Inspire 594 A Sterilization Center                     | CaviWipes™            | Compatible                  |
| Inspire 594 A Sterilization Center                     | Envirocide™           | Compatible                  |
| Inspire 594 A Sterilization Center                     | CaviCide1™            | Compatible                  |
| Inspire 594 A Sterilization Center                     | CaviWipes1™           | Compatible                  |
| Inspire 594 A Sterilization Center                     | CaviWipes™ 2.0        | Compatible                  |
| Inspire 594 A Sterilization Center                     | CaviWipes™ HP         | No information              |
| Inspire 594 A Sterilization Center                     | CaviWipes™ Bleach     | Not recommended for use     |
| Inspire 595 A Wall Mounted Cabinet/Xray Storage Insert | CaviCide™             | Compatible                  |
| Inspire 595 A Wall Mounted Cabinet/Xray Storage Insert | CaviWipes™            | Compatible                  |
| Inspire 595 A Wall Mounted Cabinet/Xray Storage Insert | Envirocide™           | Compatible                  |
| Inspire 595 A Wall Mounted Cabinet/Xray Storage Insert | CaviCide1™            | Compatible                  |

| <b>Equipment</b>                                       | <b>Metrex Product</b> | <b>Compatibility Status</b> |
|--------------------------------------------------------|-----------------------|-----------------------------|
| Inspire 595 A Wall Mounted Cabinet/Xray Storage Insert | CaviWipes1™           | Compatible                  |
| Inspire 595 A Wall Mounted Cabinet/Xray Storage Insert | CaviWipes™ 2.0        | Compatible                  |
| Inspire 595 A Wall Mounted Cabinet/Xray Storage Insert | CaviWipes™ HP         | No information              |
| Inspire 595 A Wall Mounted Cabinet/Xray Storage Insert | CaviWipes™ Bleach     | Not recommended for use     |
| A-dec 511                                              | CaviCide™             | Compatible                  |
| A-dec 511                                              | CaviWipes™            | Compatible                  |
| A-dec 511                                              | Envirocide™           | Compatible                  |
| A-dec 511                                              | CaviCide1™            | Compatible                  |
| A-dec 511                                              | CaviWipes1™           | Compatible                  |
| A-dec 511                                              | CaviWipes™ 2.0        | Compatible                  |
| A-dec 511                                              | CaviWipes™ HP         | Not recommended for use     |
| A-dec 511                                              | CaviWipes™ Bleach     | Not recommended for use     |
| A-dec Upholstery                                       | CaviCide™             | Compatible                  |
| A-dec Upholstery                                       | CaviWipes™            | Compatible                  |
| A-dec Upholstery                                       | Envirocide™           | Compatible                  |
| A-dec Upholstery                                       | CaviCide1™            | Compatible                  |
| A-dec Upholstery                                       | CaviWipes1™           | Compatible                  |
| A-dec Upholstery                                       | CaviWipes™ 2.0        | Compatible                  |
| A-dec Upholstery                                       | CaviWipes™ HP         | Not recommended for use     |
| A-dec Upholstery                                       | CaviWipes™ Bleach     | Not recommended for use     |
| Model 500 LED                                          | CaviCide™             | No information              |
| Model 500 LED                                          | CaviWipes™            | No information              |
| Model 500 LED                                          | Envirocide™           | No information              |
| Model 500 LED                                          | CaviCide1™            | No information              |

| <b>Equipment</b> | <b>Metrex Product</b> | <b>Compatibility Status</b> |
|------------------|-----------------------|-----------------------------|
| Model 500 LED    | CaviWipes1™           | No information              |
| Model 500 LED    | CaviWipes™ 2.0        | No information              |
| Model 500 LED    | CaviWipes™ HP         | No information              |
| Model 500 LED    | CaviWipes™ Bleach     | Not recommended for use     |
| Model 300 LED    | CaviCide™             | No information              |
| Model 300 LED    | CaviWipes™            | No information              |
| Model 300 LED    | Envirocide™           | No information              |
| Model 300 LED    | CaviCide1™            | No information              |
| Model 300 LED    | CaviWipes1™           | No information              |
| Model 300 LED    | CaviWipes™ 2.0        | No information              |
| Model 300 LED    | CaviWipes™ HP         | No information              |
| Model 300 LED    | CaviWipes™ Bleach     | Not recommended for use     |
| Model 421        | CaviCide™             | Compatible                  |
| Model 421        | CaviWipes™            | Compatible                  |
| Model 421        | Envirocide™           | Compatible                  |
| Model 421        | CaviCide1™            | Compatible                  |
| Model 421        | CaviWipes1™           | Compatible                  |
| Model 421        | CaviWipes™ 2.0        | Compatible                  |
| Model 421        | CaviWipes™ HP         | Not recommended for use     |
| Model 421        | CaviWipes™ Bleach     | Not recommended for use     |
| Model 422        | CaviCide™             | Compatible                  |
| Model 422        | CaviWipes™            | Compatible                  |
| Model 422        | Envirocide™           | Compatible                  |
| Model 422        | CaviCide1™            | Compatible                  |
| Model 422        | CaviWipes1™           | Compatible                  |
| Model 422        | CaviWipes™ 2.0        | Compatible                  |
| Model 422        | CaviWipes™ HP         | Not recommended for use     |

| <b>Equipment</b> | <b>Metrex Product</b> | <b>Compatibility Status</b> |
|------------------|-----------------------|-----------------------------|
| Model 422        | CaviWipes™ Bleach     | Not recommended for use     |
| Model 521        | CaviCide™             | Compatible                  |
| Model 521        | CaviWipes™            | Compatible                  |
| Model 521        | Envirocide™           | Compatible                  |
| Model 521        | CaviCide1™            | Compatible                  |
| Model 521        | CaviWipes1™           | Compatible                  |
| Model 521        | CaviWipes™ 2.0        | Compatible                  |
| Model 521        | CaviWipes™ HP         | Not recommended for use     |
| Model 521        | CaviWipes™ Bleach     | Not recommended for use     |
| Model 522        | CaviCide™             | Compatible                  |
| Model 522        | CaviWipes™            | Compatible                  |
| Model 522        | Envirocide™           | Compatible                  |
| Model 522        | CaviCide1™            | Compatible                  |
| Model 522        | CaviWipes1™           | Compatible                  |
| Model 522        | CaviWipes™ 2.0        | Compatible                  |
| Model 522        | CaviWipes™ HP         | Not recommended for use     |
| Model 522        | CaviWipes™ Bleach     | Not recommended for use     |