

Material Compatibility Status

KaVo		
Equipment	Metrex Product	Compatibility Status
All	CaviCide™	Compatible
All	CaviWipes™	Compatible
All	Envirocide™	Compatible
All	CaviCide1™	No information
All	CaviWipes1™	No information
All	CaviWipes™ 2.0	No information
All	CaviWipes™ HP	No information
All	CaviWipes™ Bleach	No information
DIAGNOCam 2170	CaviCide™	Compatible
DIAGNOCam 2170	CaviWipes™	Compatible
DIAGNOCam 2170	Envirocide™	Compatible
DIAGNOCam 2170	CaviCide1™	No information
DIAGNOCam 2170	CaviWipes1™	No information
DIAGNOCam 2170	CaviWipes™ 2.0	No information
DIAGNOCam 2170	CaviWipes™ HP	No information
DIAGNOCam 2170	CaviWipes™ Bleach	No information
DEXCam 3	CaviCide™	No information
DEXCam 3	CaviWipes™	No information
DEXCam 3	Envirocide™	No information
DEXCam 3	CaviCide1™	No information
IXS sensor & cable	CaviCide™	Compatible

Equipment	Metrex Product	Compatibility Status
IXS sensor & cable	CaviWipes™	Compatible
IXS sensor & cable	Envirocide™	Compatible
IXS sensor & cable	CaviCide1™	No information
IXS sensor & cable	CaviWipes1™	No information
IXS sensor & cable	CaviWipes™ 2.0	No information
IXS sensor & cable	CaviWipes™ HP	No information
IXS sensor & cable	CaviWipes™ Bleach	No information
PROPHYflex 4	CaviCide™	Compatible
PROPHYflex 4	CaviWipes™	Compatible
PROPHYflex 4	Envirocide™	Compatible
PROPHYflex 4	CaviCide1™	No information
PROPHYflex 4	CaviWipes1™	No information
PROPHYflex 4	CaviWipes™ 2.0	No information
PROPHYflex 4	CaviWipes™ HP	No information
PROPHYflex 4	CaviWipes™ Bleach	No information
Orthopantomograph OP 3D	CaviCide™	Compatible
Orthopantomograph OP 3D	CaviWipes™	Compatible
Orthopantomograph OP 3D	Envirocide™	Compatible
Orthopantomograph OP 3D	CaviCide1™	No information
Orthopantomograph OP 3D	CaviWipes1™	No information
Orthopantomograph OP 3D	CaviWipes™ 2.0	No information
Orthopantomograph OP 3D	CaviWipes™ HP	No information
Orthopantomograph OP 3D	CaviWipes™ Bleach	No information
DEXIS Titanium	CaviCide™	Compatible
DEXIS Titanium	CaviWipes™	Compatible
DEXIS Titanium	Envirocide™	Compatible

Equipment	Metrex Product	Compatibility Status
DEXIS Titanium	CaviCide1™	No information
DEXIS Titanium	CaviWipes1™	No information
DEXIS Titanium	CaviWipes™ 2.0	No information
DEXIS Titanium	CaviWipes™ HP	No information
DEXIS Titanium	CaviWipes™ Bleach	No information
KaVo Nomad Pro 2	CaviCide™	Compatible
KaVo Nomad Pro 2	CaviWipes™	Compatible
KaVo Nomad Pro 2	Envirocide™	Compatible
KaVo Nomad Pro 2	CaviCide1™	No information
KaVo Nomad Pro 2	CaviWipes1™	No information
KaVo Nomad Pro 2	CaviWipes™ 2.0	No information
KaVo Nomad Pro 2	CaviWipes™ HP	No information
KaVo Nomad Pro 2	CaviWipes™ Bleach	No information
Scan eXam One Plate	CaviCide™	Compatible
Scan eXam One Plate	CaviWipes™	Compatible
Scan eXam One Plate	Envirocide™	Compatible
Scan eXam One Plate	CaviCide1™	No information
Scan eXam One Plate	CaviWipes1™	No information
Scan eXam One Plate	CaviWipes™ 2.0	No information
Scan eXam One Plate	CaviWipes™ HP	No information
Scan eXam One Plate	CaviWipes™ Bleach	No information
IXS	CaviCide™	Compatible
IXS	CaviWipes™	Compatible
IXS	Envirocide™	Compatible
IXS	CaviCide1™	No information
IXS	CaviWipes1™	No information

Equipment	Metrex Product	Compatibility Status
IXS	CaviWipes™ 2.0	No information
IXS	CaviWipes™ HP	No information
IXS	CaviWipes™ Bleach	No information
Treatment Units Esthetica E70	CaviCide™	Compatible
Treatment Units Esthetica E70	CaviWipes™	Compatible
Treatment Units Esthetica E70	Envirocide™	Compatible
Treatment Units Esthetica E70	CaviCide1™	No information
Treatment Units Esthetica E70	CaviWipes1™	No information
Treatment Units Esthetica E70	CaviWipes™ 2.0	No information
Treatment Units Esthetica E70	CaviWipes™ HP	No information
Treatment Units Esthetica E70	CaviWipes™ Bleach	No information
Treatment Units Esthetica E80 Vison	CaviCide™	Compatible
Treatment Units Esthetica E80 Vison	CaviWipes™	Compatible
Treatment Units Esthetica E80 Vison	Envirocide™	Compatible
Treatment Units Esthetica E80 Vison	CaviCide1™	No information
Treatment Units Esthetica E80 Vison	CaviWipes1™	No information
Treatment Units Esthetica E80 Vison	CaviWipes™ 2.0	No information

Equipment	Metrex Product	Compatibility Status
Treatment Units Esthetica E80 Vison	CaviWipes™ HP	No information
Treatment Units Esthetica E80 Vison	CaviWipes™ Bleach	No information
Imaging all	CaviCide™	No information
Imaging all	CaviWipes™	No information
Imaging all	Envirocide™	No information
Imaging all	CaviCide1™	No information
Imaging all	CaviWipes1™	No information
Imaging all	CaviWipes™ 2.0	No information
Imaging all	CaviWipes™ HP	No information
Imaging all	CaviWipes™ Bleach	No information
i-CAT FLX V-Series	CaviCide™	Compatible
i-CAT FLX V-Series	CaviWipes™	Compatible
i-CAT FLX V-Series	Envirocide™	Compatible
i-CAT FLX V-Series	CaviCide1™	No information
i-CAT FLX V-Series	CaviWipes1™	No information
i-CAT FLX V-Series	CaviWipes™ 2.0	No information
i-CAT FLX V-Series	CaviWipes™ HP	No information
i-CAT FLX V-Series	CaviWipes™ Bleach	No information
OP 2D	CaviCide™	Compatible
OP 2D	CaviWipes™	Compatible
OP 2D	Envirocide™	Compatible
OP 2D	CaviCide1™	No information
OP 2D	CaviWipes1™	No information
OP 2D	CaviWipes™ 2.0	No information
OP 2D	CaviWipes™ HP	No information

Equipment	Metrex Product	Compatibility Status
OP 2D	CaviWipes™ Bleach	No information
DEXCam 3	CaviWipes1™	No information
DEXCam 3	CaviWipes™ 2.0	No information
DEXCam 3	CaviWipes™ HP	No information
DEXCam 3	CaviWipes™ Bleach	No information
SOREDEX DIGORA Optime	CaviCide™	Compatible
SOREDEX DIGORA Optime	CaviWipes™	Compatible
SOREDEX DIGORA Optime	Envirocide™	Compatible
SOREDEX DIGORA Optime	CaviCide1™	No information
SOREDEX DIGORA Optime	CaviWipes1™	No information
SOREDEX DIGORA Optime	CaviWipes™ 2.0	No information
SOREDEX DIGORA Optime	CaviWipes™ HP	No information
SOREDEX DIGORA Optime	CaviWipes™ Bleach	No information
DEXCam 4 HD	CaviCide™	Compatible
DEXCam 4 HD	CaviWipes™	Compatible
DEXCam 4 HD	Envirocide™	Compatible
DEXCam 4 HD	CaviCide1™	Compatible
DEXCam 4 HD	CaviWipes1™	Compatible
DEXCam 4 HD	CaviWipes™ 2.0	No information
DEXCam 4 HD	CaviWipes™ HP	No information
DEXCam 4 HD	CaviWipes™ Bleach	No information
ProXam iP	CaviWipes1™	No information
ProXam iP	CaviWipes™ 2.0	No information
ProXam iP	CaviWipes™ HP	No information
ProXam iP	CaviWipes™ Bleach	No information
ProXam iS	CaviCide™	Compatible

Equipment	Metrex Product	Compatibility Status
ProXam iS	CaviWipes™	Compatible
ProXam iS	Envirocide™	Compatible
ProXam iS	CaviCide1™	No information
ProXam iS	CaviWipes1™	No information
ProXam iS	CaviWipes™ 2.0	No information
ProXam iS	CaviWipes™ HP	No information
ProXam iS	CaviWipes™ Bleach	No information
ProXam iX	CaviCide™	No information
ProXam iX	CaviWipes™	No information
ProXam iX	Envirocide™	No information
ProXam iX	CaviCide1™	No information
ProXam iX	CaviWipes1™	No information
ProXam iX	CaviWipes™ 2.0	No information
ProXam iX	CaviWipes™ HP	No information
ProXam iX	CaviWipes™ Bleach	No information
Esthetica E50	CaviCide™	No information
Esthetica E50	CaviWipes™	No information
Esthetica E50	Envirocide™	No information
Esthetica E50	CaviCide1™	No information
Esthetica E50	CaviWipes1™	No information
Esthetica E50	CaviWipes™ 2.0	No information
Esthetica E50	CaviWipes™ HP	No information
Esthetica E50	CaviWipes™ Bleach	No information
ProXam	CaviCide™	No information
ProXam	CaviWipes™	No information
ProXam	Envirocide™	No information

Equipment	Metrex Product	Compatibility Status
ProXam	CaviCide1™	No information
ProXam	CaviWipes1™	No information
ProXam	CaviWipes™ 2.0	No information
ProXam	CaviWipes™ HP	No information
ProXam	CaviWipes™ Bleach	No information
ProXam iP	CaviCide™	No information
ProXam iP	CaviWipes™	No information
ProXam iP	Envirocide™	No information
ProXam iP	CaviCide1™	No information

SDI		
Equipment	Metrex Product	Compatibility Status
Radii Plus	CaviCide™	Not recommended for use
Radii Plus	CaviWipes™	No information
Radii Plus	Envirocide™	Not recommended for use
Radii Plus	CaviCide1™	Not recommended for use
Radii Plus	CaviWipes1™	No information
Radii Plus	CaviWipes™ 2.0	No information
Radii Plus	CaviWipes™ HP	No information
Radii Plus	CaviWipes™ Bleach	No information