

# Material Compatibility Status

MediStim		
Equipment	Metrex Product	Compatibility Status
PV	CaviCide™	Compatible
PV	CaviWipes™	Compatible
PV	Envirocide™	Compatible
PV	CaviCide1™	No information
PV	CaviWipes1™	No information
PV	CaviWipes™ 2.0	No information
PV	CaviWipes™ HP	No information
PV	CaviWipes™ Bleach	No information
PD	CaviCide™	Compatible
PD	CaviWipes™	Compatible
PD	Envirocide™	Compatible
PD	CaviCide1™	No information
PD	CaviWipes1™	No information
PD	CaviWipes™ 2.0	No information
PD	CaviWipes™ HP	No information
PD	CaviWipes™ Bleach	No information

Posey		
Equipment	Metrex Product	Compatibility Status
Bed 8060	CaviCide™	No information
Bed 8060	CaviWipes™	No information

<b>Equipment</b>	<b>Metrex Product</b>	<b>Compatibility Status</b>
Bed 8060	Envirocide™	No information
Bed 8060	CaviCide1™	No information
Bed 8060	CaviWipes1™	No information
Bed 8060	CaviWipes™ 2.0	No information
Bed 8060	CaviWipes™ HP	No information
Bed 8060	CaviWipes™ Bleach	No information
Bed 8070	CaviCide™	Compatible
Bed 8070	CaviWipes™	Compatible
Bed 8070	Envirocide™	Compatible
Bed 8070	CaviCide1™	No information
Bed 8070	CaviWipes1™	No information
Bed 8070	CaviWipes™ 2.0	Compatible
Bed 8070	CaviWipes™ HP	Compatible
Bed 8070	CaviWipes™ Bleach	Compatible
Floor cushions	CaviCide™	Compatible
Floor cushions	CaviWipes™	Compatible
Floor cushions	Envirocide™	Compatible
Floor cushions	CaviCide1™	Compatible
Floor cushions	CaviWipes1™	Compatible
Floor cushions	CaviWipes™ 2.0	Compatible
Floor cushions	CaviWipes™ HP	Compatible
Floor cushions	CaviWipes™ Bleach	Compatible

**Puritan Bennett**

<b>Equipment</b>	<b>Metrex Product</b>	<b>Compatibility Status</b>
Breeze sleep Gear	CaviCide™	No information

<b>Equipment</b>	<b>Metrex Product</b>	<b>Compatibility Status</b>
Breeze sleep Gear	CaviWipes™	No information
Breeze sleep Gear	Envirocide™	No information
Breeze sleep Gear	CaviCide1™	No information
Breeze sleep Gear	CaviWipes1™	No information
Breeze sleep Gear	CaviWipes™ 2.0	No information
Breeze sleep Gear	CaviWipes™ HP	No information
Breeze sleep Gear	CaviWipes™ Bleach	No information
Model 7200	CaviCide™	No information
Model 7200	CaviWipes™	No information
Model 7200	Envirocide™	No information
Model 7200	CaviCide1™	No information
Model 7200	CaviWipes1™	No information
Model 7200	CaviWipes™ 2.0	No information
Model 7200	CaviWipes™ HP	No information
Model 7200	CaviWipes™ Bleach	No information
Model 840	CaviCide™	No information
Model 840	CaviWipes™	No information
Model 840	Envirocide™	No information
Model 840	CaviCide1™	No information
Model 840	CaviWipes1™	No information
Model 840	CaviWipes™ 2.0	No information
Model 840	CaviWipes™ HP	No information
Model 840	CaviWipes™ Bleach	No information

### **Tektronix**

<b>Equipment</b>	<b>Metrex Product</b>	<b>Compatibility Status</b>
Model 414	CaviCide™	No information

<b>Equipment</b>	<b>Metrex Product</b>	<b>Compatibility Status</b>
Model 414	CaviWipes™	No information
Model 414	Envirocide™	No information
Model 414	CaviCide1™	No information
Model 414	CaviWipes1™	No information
Model 414	CaviWipes™ 2.0	No information
Model 414	CaviWipes™ HP	No information
Model 414	CaviWipes™ Bleach	No information